

Product datasheet for **MR212524L4V**

Fam19a3 (Tafa3) (NM_183224) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Fam19a3 (Tafa3) (NM_183224) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Tafa3 |
| Synonyms: | 7530404M11Rik; Fam19; Fam19a3; Sa; Sam3; Taf; Tafa; Tafa-3 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_183224 |
| ORF Size: | 399 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR212524). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_183224.2 , NP_899047.1 |
| RefSeq Size: | 4317 bp |
| RefSeq ORF: | 399 bp |
| Locus ID: | 329731 |
| UniProt ID: | Q7TPG6 |
| Cytogenetics: | 3 F2.2 |


[View online »](#)

Gene Summary:

This gene is a member of the Tafa family which is composed of five highly homologous genes that encode small secreted proteins. These proteins contain conserved cysteine residues at fixed positions, and are distantly related to MIP-1alpha, a member of the CC-chemokine family. The Tafa proteins are predominantly expressed in specific regions of the brain, and are postulated to function as brain-specific chemokines or neurokinins, that act as regulators of immune and nervous cells. [provided by RefSeq, Jul 2008]