

Product datasheet for **MR212616A1V**

Mouse Olfr819 (NM_001165944) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr819 (NM_001165944) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr819
Synonyms:	MOR265-2_p; Olfr819-ps1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR212616 representing NM_001165944 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCTAACTCTGAGCTCATGAAGAATGGAAGCCTGTCCCTCTGCACTGAATTCACATTGGTGGCCTTTT
CTTCTCTAGCAGAGCTGCAGCTTGCTCTTTGTTGTGTTCTTGGTTCTCTACTTGTTCACAGTGGGAGG
AAACCTCACCATCATCTGTGTGATCTGGACCACACCTTCTCTACACACGCCCATGTACTTCTTCTGGCC
AACCTCTCCTTCTAGAGATGTGCTATATCAGCAGTGTGGTGCCTCAGATGCTGGTGCACCTGCTGGTAC
AGCTCAAGACCATAAGTGTGGCAGGCTGTGCAGCTCAGATGTATGATTACCATCTTGGGGCTGACAGA
GTGTTGCCTTCTGGCAACTATGGCTTATGATCGCTTTGTGGCTATTTGTTACCCACTACATTACACTCTG
CGGATGGACCCTTCTGTCTGTCTGAAGTTGGCTGGTGCATCTTGGATGACTGGGATCCTGGTAGAATCAG
CCCAGACCACCTGGATCTTCACTCTGCCCTTTTGTGGAGCAGGAACAATTCAACACTTTTTCTGTGACAT
CATGCCTGTGGTAAAATTAGCTTGTGTGGATACTTCCAGAATGAAATTGTGATATTTATTATTTCTTTT
ATCTTCATTATGAGTCCCTGTCTCTTCTTCTATGCTCTTATGTCCGATTATTCTGACATTCTTGAAAA
TGCCCTCAGCAGCAGGCAGACACAAAGCTTTCTCCACCTGCTCCTCTCATATCCTGGTTGTTCTCTTTT
CTATGGCACTGCCTTATTCATTATCTCAACCAAGAGTTCACACACCCAGACACAGACAAGGTTACT
GCTCTCATGTACTGTGGTCAACCCTGCTCTGAATCCTGTTATCTATACGTTGAGGAACAAAGAAGTGA
AGGAAGCCTTTCAAAGGTAACACAGAGGAACTTCATAGACAAGAGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR212616 representing NM_001165944 Red=Cloning site Green=Tags(s) MSNSELMKNGSLSLCTEFTLVAFSSLAELQLVLFVVFLVLYLFTVGGNLTIIICVIWTTPSLHTPMYFFLA NLSFLEMCYISSVVPQMLVHLLVQLKTISVAGCAAQMYVFTILGLTECCLLATMAYDRFVAICYPLHYTL RMDPSVCLKLAGASWMTGILVESAQTTWIFTLPCGAGTIQHFFCDIMPVVKLACVDTSQNEIVIFIISF IFIMSPCLFILCSYVRIILTFLKMPAAGRHKAFSTCSSHILVVSIFYGTALFTYLQPKSSHTPDTDKVT ALMYTVVTPALNPVIYTLRNKEVKEAFQKVTQRKLHRQRD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001165944
ORF Size:	960 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001165944.1, NP_001159416.1</u>
RefSeq Size:	963 bp
RefSeq ORF:	963 bp
Locus ID:	100043200
Cytogenetics:	10 D3
MW:	35.9 kDa