

## Product datasheet for **MR212618**

### **Olfr361 (NM\_146368) Mouse Tagged ORF Clone**

#### Product data:

|                           |                                                                                           |
|---------------------------|-------------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                       |
| Product Name:             | Olfr361 (NM_146368) Mouse Tagged ORF Clone                                                |
| Tag:                      | Myc-DDK                                                                                   |
| Symbol:                   | Olfr361                                                                                   |
| Synonyms:                 | MOR159-3                                                                                  |
| Mammalian Cell Selection: | Neomycin                                                                                  |
| Vector:                   | pCMV6-Entry (PS100001)                                                                    |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                                      |
| ORF Nucleotide Sequence:  | <p>&gt;MR212618 representing NM_146368</p> <p>Red=Cloning site Blue=ORF Green=Tags(s)</p> |

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCTACTAAAAACAGGACAGAAAGTGACTGAATTTGTGTTGTTGGGATTGTCCAGCCAGCCAGAGATGC  
 AGCCAGTTATTTTGGAGTTGTCCTTATCATGTACCTGATGGCAGTTTGGGCAATACCCTACTAGTGCT  
 TGTTGCTTGCTCAGACCCCAAGCTTCATACCCCATGACTTCCTTCTCAGTCAGCTTTCCTTGATTGAC  
 ATATCTCTAACTATCATTGTTCCCGAGATGCTGGTGACACACTTCTGTAACAGAACCATTTCT  
 ATAAGTCTGCATGACACAGCTCTTTCTTTATGACAGTGGGAGTATGGAAGTTCAGTTGCTTGGTGC  
 AATGGCATATGACCGCTATGTTGCCATCTGTGACCCCTTAAGATATTCTCCATTGTGAGCCATAGCCTC  
 TGTCTGCGAATAACATTGACCTCATGGGTGGTGGTCAGCCTTAACAGTCTCCTGTACAGTGTGTTGGTCA  
 CTCGTTTAACTTCTGTGGCAACAAGGTCACCCATTTCTTTGTGACATCACTCCTCTACTGAAGCTCTC  
 CTGTACTCGACCAAGTGGTCAATGAAATGTTAATCTTCACTGAGGGTGTAGCTGTGGTGGTAAGCCCCATT  
 TTCTTTATTTTGGCTCCTATGTTGCGATAGGTATTGTTATAGCCCATGCACTCAGTTGCTGCCCTGC  
 AGAAAGCTCTGTCCACTGTAGCTCTCATGTTCTTGTGTGATGTTTCTGTTGGTTCTTTGGTTACAT  
 GTACCTCAAGCCATCATCCAGCTACAATCTAGAACAAGATCGCCAGGTTGCCATCTTTACACACTCATT  
 TCTCCCATGTTAAATCCACTAATCTACAGCCTCAGGAACAGGAAGTTAAAGGAGCTCATGGAGGCTTT  
 TTAGGAAACTATATCTCAGGCAGTGTCCAGCTTGGATTCCCAACACATGCCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

**Protein Sequence:** >MR212618 representing NM\_146368  
 Red=Cloning site Green=Tags(s)

MATKNRTEVFVLLGLSSQPEMQPVIFGVVLIYMLMAVLGNTLLVLVACSDPKLHTPMYFLLSQLSLID  
 ISLTTIIVPQMLVHTLSVNRITISYNCMTQLFSFMTVGSMEVHLLGAMAYDRYVAICDPLRYSSIVSHSL  
 CLRITLTSWVVVSLNLLYSVLVTRLFCGNKVTHFFCDITPLKLCTRPVVMNEMIFTEGVAVVVSPI  
 FFIFGSYVRIGIVIAHMSVAALQKALSTCSSHVLVVMFLFGSLVHMYLKPSSSYNLEQDRQVAIFYTLI  
 SPMLNPLIYSLRNQEVKGALWRLFRKLYLRQCPAWISHNTCQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9083\\_f12.zip](https://cdn.origene.com/chromatograms/mm9083_f12.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_146368

**ORF Size:** 966 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_146368.1](#), [NP\\_666480.1](#)

**RefSeq Size:** 969 bp

**RefSeq ORF:** 969 bp


**Locus ID:** 258365

**Cytogenetics:** 2 B

**MW:** 36.2 kDa

**Gene Summary:** Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

## Product images:



Circular map for MR212618