

Product datasheet for MR212619

Ceacam19 (NM_177036) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ceacam19 (NM_177036) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ceacam19
Synonyms: C130022P09Rik; CEAL1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR212619 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAGACCTGGACCCACCGGGCTGCTTCTCAAAGGGCCTCCTGTTCTCAGCCTTGTTCTGGCCCTTG
AGCTGCCCAAGGCTCACAGGCAGCCCTCCACATCGAGATGATCCAGAGCAGCCTCAAAGAACCAGGA
CCTTCTTCTGTCCCTTACGGTGTCCAGCACTGCTCAGGATTCATCTGGTACCTGGGGGAGGAGGCT
TCTGGAGGCACAAGGCTGTTCTTATATCCCCGACTACCACGGCCCCAGAGGGATGGCAGTGCCATGG
GACAGCGAGACATTGTTGGCTTCCAAAATGGTTCCATGCTGCTTCGAAATGCCAACCTTCAGATAGTGG
CACCTATCAGGTGGCCGCCACCATGAACCTGCCTGGACAATGAAGGCCAAGATTAACGTCCATGTAGCT
GGTGCACGCGAGGTCTCCTCACTCCTCCCTGTCCATGCCGGGATCACGGTAGCCATCATATTGGGAGTC
TTGCTATCGGGTCACTACTTGTATGTGGCATTGCCTATGTCTGGTCACTCGAAGCCGGAAGGGCCGGAG
CCCCAGGAATCCAGCCACAGAGAAGCCAGAAGTGAACCTGAGTCTGAATCTGGTGACAGCAACATCTAT
GAAGTGATGCCATCTCCAGTCTTCTGGTGTCTCCATCAGTGACATGGAGCAAGGAAGCCACCCATGA
GTCAGCCTCCACCACCTCCACCAGCATCTCAGGAGCCAGAGCCAGAGCACCAGCACTATCAGGACCTACT
AAACCTGACCCAGCCCCATACTGCCAGCTGGTGCCAACCTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_177036.5](#), [NP_796010.1](#)

RefSeq Size: 3361 bp


RefSeq ORF: 816 bp

Locus ID: 319930

Cytogenetics: 7 A3

MW: 29.1 kDa

Product images:



Circular map for MR212619