

Product datasheet for MR212709

AA792892 (NM_178894) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AA792892 (NM_178894) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	AA792892
Synonyms:	MGC159137
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR212709 representing NM_178894 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGTGGTCAGACCCACCCACACGCACTCTGCTGACACATGAAGCTTTGATCAAGTCTCTCTGGAGG
AGTTGCCACTGTGGTCTTCCAGCACTATTCAAAGAGGTCTTTGCTGGCAGACACACCAATCTCTAAA
GGCAATGGTGGCAGCCTGGCCTTCCCCTGTCCCTGTGGGACCATTGATAAAGAAGCCTAACCTAGAA
ACCTTGAAGCTCTGCTAGATGGAGTAGACATGAGACTGACAAGAGATTTACCCAGGTGCCTGATGA
CTCCCTTCGAGACCTTCCATCACCAATGCCTAATTTACAGAGAGACTTTGATTTCTTTCTTGCTC
CCACAACCTTTTAAATTAACATCTGAAATAAGAGGGATGATCTTATATGCTTTGGATCTTATGCC
CTGAGAGGTCTCTACAAAAAGTAGCAGACACTCTTGAGTTTCTGGATTTTCAGTGGTGTAGGATGAAGG
ACTCCAGATCAATGCCCTCCTACCTGCCCTCAGTCAATGCACTCACCTCAACCAGATCAACTTCTACAA
CAATGACTTCTCCATGTCCACCCGGAAGGAGCTTTACAGCACACAGCCAAGTGGAGCAAGATGAATGTG
GAACAATGCCCTGCCCTCTGCAGTGTATAATGAGTTGGGTGATGTCTCTGTAGAAAGATTTGTCCAAC
TTTGTGAGGAGCTCATGTACTCTCAGAGCTATAAGGCAGCCCAAGAGCATCTTTTGTACAGATAG
CTGTACAAAATGTTGTAAGCCCTGTGTCTATGACCAGGACCCTATACTTTGTTCTTGTTTTCAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >MR212709 representing NM_178894
Red=Cloning site Green=Tags(s)

MSGQTPPTRLLTHEALIKSSLEELPTVVPALFKEVFAGRHTNLLKAMVAAWPFCLPVGPLIKKPNLE
 TLQALLDGVDMRLTREFHPRCLMTPFETFSITQCLISQRDFDFSCSHNLFKLEIRGMILYALDLMP
 LRGLLQKVADTLEFLDFQWCRMKDSQINALLPALSQCTHLNQINFYNDFSMSTLKELLQHTANWSKMNV
 EQCPAPLQCYNELGHVSVRFVQLCQELMYTLRAIRQPKSISFATDSCHKCKPCVYDQDPILCSCFQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9062_f08.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_178894

ORF Size: 834 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_178894.4](#), [NP_849225.2](#)

RefSeq Size: 1212 bp


RefSeq ORF: 837 bp

Locus ID: 100554

Cytogenetics: 5 E3

MW: 31.9 kDa

Product images:



Circular map for MR212709