

Product datasheet for **MR212709A1V**

Mouse AA792892 (NM_178894) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse AA792892 (NM_178894) AAV Particle
Tag: Myc-DDK
Symbol: AA792892
Synonyms: MGC159137
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR212709 representing NM_178894
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAGTGGTCAGACCCACCCACACGCACTCTGCTGACACATGAAGCTTTGATCAAGTCCTCTCTGGAGG
 AGTTGCCTACTGTGGTCTTCCCAGCACTATTCAAAGAGGTCTTTGCTGGCAGACACCAATCTCCTAAA
 GGCAATGGTGGCAGCCTGGCCTTTCCCCTGTCTCCCTGTGGGACCATTGATAAAGAAGCCTAACCTAGAA
 ACCTTGCAAGCTCTGCTAGATGGAGTAGACATGAGACTGACAAGAGAGTTTACCCCAGGTGCCTGATGA
 CTCCCTTCGAGACCTTCTCCATCACCAATGCCTAATTCACAGAGAGACTTTGATTTCTTTCTTGCTC
 CCACAACCTCTTTAAATTAACATCTGGAATAAGAGGGATGATCTTATATGCTTTGGATCTTATGCCT
 CTGAGAGGTCTCTACAAAAAGTAGCAGACACTCTTGAGTTTCTGGATTTTCAGTGGTGTAGGATGAAGG
 ACTCCCAGATCAATGCCCTCTACCTGCCCTCAGTCAATGCACTCACCTCAACCAGATCAACTTCTACAA
 CAATGACTTCTCCATGTCCACCCTGAAGGAGCTTTACAGCACACAGCCAAGTGGAGCAAGATGAATGTG
 GAACAATGCCCTGCCCTCTGCAGTGTATAATGAGTTGGTTCATGTCTCTGTAGAAAGATTTGTCCAAC
 TTTGTGAGGAGCTCATGTATACTCTCAGAGCTATAAGGCAGCCCAAGAGCATCTTTTGTCTACAGATAG
 CTGTACAAATGTTGAAGCCCTGTGTCTATGACCAGGACCCTATACTTTGTTCTGTTTTCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>MR212709 representing NM_178894 Red=Cloning site Green=Tags(s) MSGQTPPTRTLLTHEALIKSSLEELPTVVFPALFKEVFAGRHTNLLKAMVAAWPFCLPVGPLIKKPNLE TLQALLDGVDMRLTREFHPRCLMTPFETFSITQCLISQRDFDFSCSHNLFKLKHLEIRGMILYALDLMP LRGLLQKVADTLEFLDFQWCRMKDSQINALLPALSQCTHLNQINFYNDFSMSTLKELLQHTANWSKMNV EQCPAPLQCYNELGHVSVERFVQLCQELMYTLRAIRQPKSISFATDSCHKCKPCVYDQDPILCSCFQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_178894
ORF Size:	834 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_178894.4, NP_849225.2</u>
RefSeq Size:	1212 bp
RefSeq ORF:	837 bp
Locus ID:	100554
Cytogenetics:	5 E3
MW:	31.9 kDa