

Product datasheet for MR212767

Olfr933 (NM_146441) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr933 (NM_146441) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr933
Synonyms:	MOR171-22
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR212767 representing NM_146441 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGGCACTGGGAATCATTCTGTAAACAGTTGTGTTTGTCTTGGTGGGATTAACACAGCAGCCTGAACTCC
TGCTGCCCTCTTCATTCTGTTCTGGGCATCTATGTGGTGACAGCAGTGGGAACCTGGGCATGATCCT
GCTCATCACAGTCAGTCCACTGCTGCATACTCCTATGTACTATTCCTCAGCAGCTTATCCTGTGTTGAT
CTCTGCTATTCCACTGTCATTACACCCAAATGCTGGTGAACCTTCTTGGGAAGAAAAATTTGATTGTCT
ACTCAGAGTGCATGGCCAGCTCTTTTCTTTGTGATCTTCGTGGTGGCTGAGGGCTACCTGCTGACTGC
CATGGCATATGATCGCTATGTGGCCATCTGTAGACCATTGCTGTATAATGTGATCATGTCCTCCAGACTT
TGCTCACTGCTAGTGTTGGTTGCTTTCATCTTAGGCTTTGTTTCTGCCTTAGCACACACAAGTGCTATGA
TGAATCTGAGCTTCTGTAAATCCACGTCATAAGCCATTATTTCTGTGATGTTCTTCCCCTCCTAAACCT
CTCCTGCTCTGACATAAAGCTGAATGAGCTTCTTCTATTCATCATTGCAGGCTTCAACACCTTGGTGCC
ACTCTAGCTGTTGCCATCTTTATGTCTTTATTTCTGCAGCATCCTTCACATCAAGTCATCAAAGGGCA
GGTCCAAAGCCTTTGGAACCTGCAGCTCTCATCTCATGGCTGTGGGGATCTTCTTTGGGTCTATCACCTT
CATGTATTTCAAGCCTCCTTCAAGTAACTCTCTGGAGCAAGAGAAGGTGTCTTCTGTGTTCTATACCACA
GTGATCCCATGCTGAACCCATTAATATACAGTTTGAGGAACAAAGATGTGAAGAAAGCATTGGGCAAAAT
GTCTGGCAGGGAGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR212767 representing NM_146441
 Red=Cloning site Green=Tags(s)


MGTGNHSTVTVFVLVGLTQQPELLLPFILFLGIYVVTAVGNLGMILLITVSPLLHTPMYYFLSSLSCVD
 LCYSTVITPKMLVNFLGKKNLIYSECAQLFFFVIFVVAEGYLLTAMAYDRYVAICRPLLNVIMSSRL
 CSLLVLVAFILGFVSALAHTSAMNLSFCKSHVISHYFCDVLPDLLNLSCSDIKLNELLLFIAGFNTLVP
 TLAVAISYVVFIFCSILHIKSSKGRSKAFGTCSSHLMVAGIFFGSITFMFKPPSSNSLEQKVSSVFYTT
 VIPMLNPLIYSLRNKDVKKALGKCLAGR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9083_h05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_146441

ORF Size: 924 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_146441.1](#), [NP_666652.1](#)

RefSeq Size: 927 bp

RefSeq ORF: 927 bp


Locus ID: 258433

Cytogenetics: 9 A5.1

MW: 34 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR212767