

Product datasheet for **MR212767A1V**

Mouse Olfr933 (NM_146441) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr933 (NM_146441) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr933
Synonyms:	MOR171-22
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR212767 representing NM_146441 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGGCACTGGGAATCATTCTGTAACAGTTGTGTTTGTCTTGGTGGGATTAACACAGCAGCCTGAACTCC
TGCTGCCCCCTTTCATTCTGTTCTGGGCATCTATGTGGTGACAGCAGTGGGGAACCTGGGCATGATCCT
GCTCATCACAGTCAGTCCACTGCTGCATACTCCTATGTACTATTTCTCAGCAGCTTATCCTGTGTTGAT
CTCTGCTATTCCACTGTCATTACACCCAAAATGCTGGTGAACCTTCTTGGGAAGAAAAATTTGATTGTCT
ACTCAGAGTGCATGGCCAGCTCTTTTCTTTGTGATCTTCGTGGTGGCTGAGGGCTACCTGCTGACTGC
CATGGCATATGATCGCTATGTGGCCATCTGTAGACCATTGCTGTATAATGTGATCATGTCCTCCAGACTT
TGCTCACTGCTAGTGTTGGTTGCTTTTCATCTTAGGCTTTGTTTCTGCCTAGCACACACAAGTGCTATGA
TGAATCTGAGCTTCTGTAATCCACGTCATAAGCCATTATTTCTGTGATGTTCTTCCCCTCCTAAACCT
CTCCTGCTCTGACATAAAGCTGAATGAGCTTCTTATTATCATCATTCAGGCTTCAACACCTTGGTGCCC
ACTCTAGCTGTTGCCATCTCTTATGTCTTTATTTCTGCAGCATCCTTCAATCAAGTCATCAAAGGGCA
GGTCAAAGCCTTTGGAACCTGCAGCTCTCATCTCATGGCTGTGGGATCTTCTTTGGGTCTATCACCTT
CATGTATTTCAAGCCTCCTTCAAGTAATCTCTGGAGCAAGAGAAGGTGTCTTCTGTGTTCTATACCACA
GTGATCCCCATGCTGAACCCATTAATATACAGTTTGAGGAACAAAGATGTGAAGAAAGCATTGGGCAAT
GTCTGGCAGGGAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR212767 representing NM_146441 Red=Cloning site Green=Tags(s) MGTGNHSVTVFVLVGLTQQPELLLPFILFLGIYVVTAVGNLGMILLITVSPLLHTPMYYFLSSLSCVD LCYSTVITPKMLVNFGLGKKNLIYSECMAQLFFVIFVVAEGYLLTAMAYDRYVAICRPLLYNVIMSSRL CSLLVLVAFILGFVSALAHTSAMMNL SFCKSHVISHYFCDVLP LLNLSCSDIKLNELLFIAGFNTLVP TLAVAI SYVFI FCSILHIKSSKGRSKAFGTCSHLM AVGIFFGSITFMYFKPPSSNSLEQEKVSSVFYTT VIPMLNPLIYSLRNKDVKKALGKCLAGR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_146441
ORF Size:	924 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_146441.1, NP_666652.1</u>
RefSeq Size:	927 bp
RefSeq ORF:	927 bp
Locus ID:	258433
Cytogenetics:	9 A5.1
MW:	34 kDa