

Product datasheet for **MR212841**

Fam169a (NM_001146045) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fam169a (NM_001146045) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Fam169a
Synonyms:	B230112C05Rik; SLAP75
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR212841 representing NM_001146045
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCATTCCCTGTGGATTTGCTGGATAACTGCACGCATGAGGAATTGGAAAACCTTTCTGAGGATTACC
TCTCCAGCCTAAGGTGTGGGGACCCTGAACATCCAGAGTGCTTTTCTAGTCTCAACATTACGATTCCCGT
TAGTCTGTCAAATGTGGGCTTTGTCCTCTCTATGGCGGCAATCAGACACAGAAAATTTTGCTCTCTTT
GCGCCTGAAGATTCTCTTACAGCTGTGGCACTTTACCTTGTGGTGTGAGTGGGCCATTGATGATTTG
TGAAAACCTCTGAACCTTCTAGAGAGGGCTGAAGCAGGTGAGCACTCTTGGAGAGAGAGTGGTTCTGTA
TGTTCTCAATCGGATTATTTATAGAAAGCAGGAAATGGAGAGAAATGAGATCCCATTCTGTGCCATAGC
AGTACTGACTATGCTAAGATTCTGTGGAAGAAAGGAGAGGCCGTTGGGTTTTATTAGTTAAGCCTACAG
GAAGCTTATGTGCCTCATTCTTACTCAGAATTACCAACTGCCAGTTCTTGATACAATGTTTATAAGGAA
GAAATACAGAGGAAAAGATTTGGGGCTTCATATGCTGGAAGATTTGTTGATTCTTTACAGAAGATGCA
CTTGGCTTACGGTACCCACTGTCTTCTCATGTACACAGCTTCCAAGCAATACTTTGAAAAGTATCCAG
GAGATCATGAACCTTTTGGGAAGTTGAAGGTGTGGGTCTTGGCATCAGCGAGTACCAGTACCCCGAGC
ATTGCAGAGAGAGGCTATTAAGCTACAGATGTTTCTCAATATGAAGCCACAAGACTGTGTCTGGAGAA
TATGGTCTTGCTGCTGTACCAGAACATGAACCTGGACTGGACGACACTCAGTCTAGTGAGCTGCAGATCC
ATTCCTAAAAGATGCCTTTGCTAGCACTTCTGAAGGTCCAGAAAAGACTCCTGTTTCCACGCGTACTCG
GAGCAGCCATCTAAAGCGGCCAAAGATTGGGAAGCATTTCAGGATTCTGAGTTTAGTATTCTCAAGGT
GAAGATGAAAACGTTGCCAAGACTTCACCTACGGCTTCAGTAAACAAAATAGAGTATGCTGCACGAAAT
CAGAGAGCTCAGAAGAATTCCTTGAAGAAGAGCCTGAGCAGGGAGTGATTGACTTTGAGGATGAAAGCGG
TGATAAAGATGCACAGCCTGCTCTGAAAACCTCAGCCACGCTGCAGAAACAAGATGGGGACAAGGACTCT
CGGTTAGAGCCTGTGAATGGGGAGGTAATGGATGCTGCCCTCAAGCCCTCTTTACAACCTGAAGACGAGG
ACTCCACCAGTGAGGGTCTTGAAGAAGACTTAAAGGTTCCACCTTTTAATTCCAGCGGAGAACCTGGAAA
TCCTGTTCCACTAGTAGCAGAGTCTTCAAAGTGCCCGAGGCCACATTAGCTAAGACCTCACCTGATACT
GACTCAGAAATGTTGATCGATCAAAGCCCATCTGATGATAAAGGACATACTGAAGAAAATCTGTCCCAG
TTTCAAAAAGAAAACACTTTTGGGGAGCTCAGACAATGTTGCCACTGTATCAAATATAGAAAAGTCTGA
TGGAATTTTCCAATTCTGTGGTACCAGAATTTCCCGAAGAACCAGTCTCTCAGAATTTATCACCCAAC
ACTACTTCTCAGTGAAGACCAGGGTGTAGGAGGAGCACCTGAGGCCAGGAACCATCGGCTACTCAGA
GTTCTCTGATAGAGGTAGAACTTGAAGACGCGCCATTTCCACAGAATGCAGGACAGAAGGCCAGTCAGA
GGAACAGTCTGAAGCGTCTCTGAGCACCTGGAACAGTTTACACAGTCAGCAGAGAAAGGTGTGGATAGC
AGCTCAGAGGAGATTGAAGTGAAGTTCCTGTAGTAGACAGGCGCAATCTGAGAAGAAAGGCCAAAGGGC
ATAAGGGCCCTGGCAAGAAGAAAGCCAAGCTGACC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR212841 representing NM_001146045
 Red=Cloning site Green=Tags(s)


MAFPVDLLDNCTHEELENSSELYLSSLRCDPEHPECFSSLNITIPVSLSNVGFVPLYGGNQTKILALF
 APEDSLTAVALYLVGQWVAIDDIIVKTSEPSREGLKQVSTLGERVLYVLNRIIYRKQEMERNEIPFLCHS
 STDYAKILWKKGEAVGFYSVKPTGSLCASFLTQNYQLPVLDTMFIRKKYRGKDLGLHMLEDFVDSFTEDA
 LGLRYPLSSLMYTASKQYFEKYPGDHELLWEVEGVGHWHQRPVTRALQREAIKATDVSQYEATRPVSGE
 YGLAAVPEHEPGLDDTQSSELIHSLKDAFASTSEGPEKTPVSTRTRSSHLLKRPKIGKHFQDSEFSSSQG
 EDENVAKTSPTASVINKIEYAARTSESSEEFLEEEPEQGVIDFEDES GDKDAQPALETQPRLQKQDGDKDS
 ALEPVNGEVMDAALKPSSLTTEDEDSTSEGLEEDLKVPPFNSSGEPGNPVPLVAESSKVPEATLAKTSPDT
 DSEMLIDQSPSDDKGHTTEENLSPVSKKTLGSSDNVATVSNIEKSDGNFPNSVVPFPEEPVSNLSPN
 TTSSVEDQGEEGAPEAQEPSATQSSLIEVELEDAPFPQAGQKSQSEEQSEASSEHLEQFTQSAEKAVDS
 SSEEIEVEVPVDRRNLRRKAKGHKGP GK K K K A K L T

SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9045_a05.zip

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_001146045

ORF Size: 1995 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146045.1](#), [NP_001139517.1](#)

RefSeq Size: 4537 bp

RefSeq ORF: 1998 bp


Locus ID: 320557

UniProt ID: [Q5XG69](#)

Cytogenetics: 13 D1

MW: 73.7 kDa

Product images:



Circular map for MR212841