

Product datasheet for **MR212890A1V**

Mouse Olfr90 (NM_146477) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr90 (NM_146477) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr90
Synonyms:	bM573K1.2; MOR256-21
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR212890 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTCAACCAGAGCTCCCCGTGGGCTTCTCTGCTGGGCTTCTCTGAACACCCACAACCTGGAAAAGG
TTCTCTTCGTCGTTGTCTTGCTCCTACCTCCTCACCCTCCTAGGAAACACACTAATCCTCCTGCTGTC
CACCTGGATCCCAGGCTCCACTCTCCAATGTACTTCTTCTCTTAACCTCTCCTTCTTGACCTCTGC
TTCACCACAACCTGTGTCCCCAGATGCTGTCAACCTCTGGGGCCCCGAAAAGACCATCAGCTTCCTGG
GATGCTTTGTCCAGCTTTCATCTTCTTATCTCTGGGGACGACGAGTGCATCCTGCTAGCAGTGATGGC
CTTTGACCGCTATGTGGCTGTGTGCCAGCCCTGCACTATGCCACCATCATCCACCTCGCCTGTGTTGC
CAGCTGGCAGCTGTGGCCTGGACGATTGGCCTAGTGGAGTCAGTAGTTCAGACGCCATCCACTCTGCGCC
TGCCCTTCTGTCCCCACCATCAGGTTGACGACTTTGTGTGTGAGGTCCCTGCTTTGATTGACTGTCTGT
TGGGGACACCACCTATAATGAGATACAAATGGCTATTGCCAGTGTCTTCATCTTGGTTGTACCTCTGAGC
CTCATCCTTGTCTTTATGGTGCCATTGCCAGGCGAGTCTGAGGATAAGCTCTGCAAAGGGGCGCAGGA
AAGCTTTTGGGACCTGCTTCCACCTCATTGTGGTCACTCTTCTACAGCTCAGTCATCGTGTTTA
TCTGCAGCCCAAGAATCCTTATGCTCGAGAGAGGGCAAGTTCTTTGGTCTTCTATGCTGTGGAACT
CCTTCTCTTAACCTCTCATATATACCCTAAGGAACAAGGAAGTAAAGAGGGCATTGAGGAGGCTGCTCT
GGAAGGAAGTGAAGCCCGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR212890 protein sequence Red=Cloning site Green=Tags(s) MVNQSSPVGFLLLGFSEHPQLEKVL FVVVLC SYLL TLLGNTL ILLLSTLDPRLHSPMYFFLSNLSFLDLC FTTTCVPQMLFNLWGPEKTISFLGCFVQLFIFLSLGTTECILLAVMAFD RYVAVCQPLHYATIIHPRLCC QLA A VAWTIGLVESVQTPSTLRLPFCPHHQVDDFVCEVPALIRLSCGDTTYNEIQMAIASVFILVVPLS LILVSYGA IARAVLR ISSAKGRRKAFGTCSSHLIVVTLFYSSVIAVYLQPKNPYARER GKF FGLFYAVGT PSLNPLIYTLRNKEVKRAFRRLLWKEVKPG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_146477
ORF Size:	930 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_146477.2, NP_666688.2</u>
RefSeq Size:	1046 bp
RefSeq ORF:	933 bp
Locus ID:	258469
Cytogenetics:	17 B1
MW:	34.7 kDa