

Product datasheet for **MR212904A1V**

Mouse Olfr890 (NM_146481) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr890 (NM_146481) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr890
Synonyms:	MOR162-3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR212904 representing NM_146481 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACATTGGAAGAATGGCCTTCAGCAATGACTCTTCTGTGAAGGAGTTATCCTGCTGGGCTTAACAC
AGCAGCCAGAGCTCCAGATGCCTCTTTCTTCTCTTCTGGGAATCTATGTGGTCTCCATGGTGGGAA
CCTGGGATTGATTGTTCTGATTGTTTGAATCCTCACCTGCACACCCCATGACTACTTCTCTTCAAC
CTTCTTTTATTGATTCTGCTACTCTTCTGTCAACACCCAAATGCTGGTGGGTTTGTGAAGCAAA
ACATCATCTCTCATGCAGAGTGTATGACTCAGCTCTTTTCTTCTCTTCTTTGTTATTGATGAGTGCTA
CATTTTGACAGCAATGGCCTATGACAGATATGCTGCAATCTGTAAGCCCTGCTTTACCGGTCACCATG
TCCTATCAGGTCTGCCTCTTGATGACAGGGGGATGTATGTGATGGGGCTTGTTGGTGCCATAGCCACA
TTGTTTGCATGCTGAGACTCACCTTCTGTGAAGGCTACATCATTAAATCATTACATGTGTGACATACCC
TCTCCTGAAGCTGTCTGCACAAGCACCTACATCAATGAGCTGGTAGTTTTCATTGTTGTGGGTGCAAT
GTCATAGTTCCACATTGACTATTTTATTTCTTACACCTTGATCCTTTCAACATCCTCAGCATCCATT
CTGCAGAAGGTAGGTCAAAGCATTGAGACCTGTGGTCCCATTTCATAGCTGTTTCTTTTCTATGG
AGCTTCAGCATTATGTATCTTAAACCTTCTAGTGCATCAGTGGATGATGATAAAATATCTACCATTTT
TATACCATTGTGGGCCAATGTTGAATCCGTTTATTTACAGTTTAAGGAATAAGGATGTCCACATTGCAC
TGAGAAAACTTTGAAGAAAAGTACGTTTATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR212904 representing NM_146481 Red=Cloning site Green=Tags(s) MTLERMAFSNDSSVKEFILLGLTQQPELQMPLFFLFLGIYVSMVGNLGLIVLIVLNPHLHTPMYYFLFN LSFIDFCYSSVITPKMLVGFVKQNIISHAECMTQLFFSFFVIDECYILTAMAYDRYAAICKPLLYQVTM SYQVCLLMTGGMYVMGLVGAIHIVCMLRLTFCEGYIINHMYCDIPLLKLSCTSTYINELVVFIVVGVN VIVPTLTIFISYTLILSNILSIHSAEGRSKAFRTCGSHFIAVSLFYGASAFMYLKPSSASVDDDKISTIF YTIVGPMLNPFYISLRNKDVHIALRKLKSTFI TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_146481
ORF Size:	942 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_146481.2, NP_666692.2</u>
RefSeq Size:	945 bp
RefSeq ORF:	945 bp
Locus ID:	258474
Cytogenetics:	9 A5.1
MW:	35.5 kDa