

Product datasheet for MR212960

Olfr153 (NM_206823) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr153 (NM_206823) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr153
Synonyms:	GA_x6K02T2Q125-49033418-49034341; MOR177-5; Olfr4-2; V5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR212960 representing NM_206823 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAATATAGGAACTGTTCTTGAATGAATTCATTTTTGTGGGAGTTACCAATAACCCAGAGATGAAAG
GGACGCTGTTTACCATATTTCTCCTCATCTATCTCATCAATCTTCTGGGAAATATTGGAATGATCATTTT
GATCAGAATGGATCCCCAGCTCCACACACCCATGACTTTTTCTGAGCCACTTGCTTTCTGTGATCTT
TGCTATCCACAGCAATTGGTCCAAGATGCTGGTAGATATGTTGGCAAGAACAAGTCAATTCCTTTCT
GGGCTGTGCTGTGCAGTTCTTCAATTTCTGTACCTTTGCAGATTCTGAGTGTGTGCTGCTGGCAGTCAT
GGCCTTTGATAGATACCGGCCATTAGCAACCCCTTGCTCTACACAGCAAACATGTCAGCAGGAAGTGT
TTCATGTTTCATGGCTGGCGTTTACCTGGTGGGAACATCAGATGCTTTGATACATACAACTTTGGCATTCC
GTTTGTGTTTCTGCGGGTCTAATGAGATCAACCACTTCTTCTGTGATCTACCTCCAATTTACCTTCTCTC
CTGCTCAGATACTCAAGTAAATTACCTGGCACTATTCATATTTATGGCTTCTTGAAGTGAAGCACCATC
TCAGGAGTCTTGTCTTACTGTTACATCATCTCATCAGTCTTAAAAATCCGTTCAACTGAAGGGAGGT
TCAAAGCTTTCTCACCTGCACCTCTCACTTGACTGCAGTTGCTATTTTTCAAGAACTCTACTCTTCCAC
GTATTTTCGACCAAGTTCTTCTTATTCCCTTGATCAAGACAAAATGACATCTCTATTTTATACCTTTGTG
ATTCCTATGCTGAACCCACTGATTTATAGTTTAAGGAACAAGGATGTGAAAGAGGCCCTTGGGGAAAATGA
AAAGGAAAAGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR212960 representing NM_206823
Red=Cloning site Green=Tags(s)

MNIGNCSLNEFIFVGVTTNPEMKGTLFTIFLLIYLINLLGNIGMIILIRMDPQLHTPMYFFLSHLSFCDL
 CYSTAIGPKMLVDMFGKNKSIPFWGCAVQFFISCTFADSECVLLAVMAFDYQAISNPLLYANMSSRKC
 FMFMAGVYLVGTSDALIHTTLAFRLCFCGSNEINHFFCDLPPIYLLSCSDTQVNYLALFTIYGFLSTI
 SGVLVSYCYIISSVLKIRSTEGRFKAFSTCTSHLTAVAFQGTLLFTYFRPSSSYSLDQDKMTSLFYTLV
 IPMLNPLIYSLRNKDKVEALGKMKRKR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9080_e09.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_206823

ORF Size: 921 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_206823.1](#), [NP_996559.1](#)

RefSeq Size: 924 bp

RefSeq ORF: 924 bp


Locus ID: 110511

Cytogenetics: 2 D

MW: 34.7 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR212960