

Product datasheet for **MR212960A1V**

Mouse Olfr153 (NM_206823) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr153 (NM_206823) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr153
Synonyms:	GA_x6K02T2Q125-49033418-49034341; MOR177-5; Olfr4-2; V5
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR212960 representing NM_206823 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAATATAGGAACTGTTCTTGAATGAATTCATTTTTGTGGGAGTTACCAATAACCCAGAGATGAAAG
GGACGCTGTTTACCATATTTCTCCTCATCTATCTCATCAATCTTCTGGGAAATATTGGAATGATCATTTT
GATCAGAATGGATCCCCAGCTCCACACACCCATGTACTTTTCTGAGCCACTTGCTTTCTGTGATCTT
TGCTATTCCACAGCAATTGGTCCCAAGATGCTGGTAGATATGTTGGCAAGAACAAGTCAATTCCTTTCT
GGGGCTGTGCTGTGCAGTTCCTCATTTCTGTACCTTTGCAGATTCTGAGTGTGTGCTGCTGGCAGTCAT
GGCCTTTGATAGATACCAGGCCATTAGCAACCCCTTGCTCTACACAGCAAACATGTCCAGCAGGAAGTGT
TTCATGTTTCATGGCTGGCGTTTACCTGGTGGGAACATCAGATGCTTTGATACATACAACCTTTGGCATTCC
GTTTGTGTTTCTGCGGGTCTAATGAGATCAACCACTTCTTCTGTGATCTACCTCCAATTTACCTTCTCTC
CTGCTCAGATACTCAAGTAAATTACCTGGCACTATTCACTATTTATGGCTTCTTGAACCTGAGCACCATC
TCAGGAGTCTTGTCTTACTGTTACATCATCTCATCAGTCTTAAAAATCCGTTCAACTGAAGGGAGGT
TCAAAGCTTTCTCCACCTGCACCTCTCACTTGACTGCAGTTGCTATTTTCAAGGAACTCTACTCTTCAC
GTATTTTCGACCAAGTCTTCTTATTCCCTTGATCAAGACAAAATGACATCTCTATTTTATACCCTTG
ATTCCTATGCTGAACCCACTGATTTATAGTTTAAGGAACAAGGATGTGAAAGAGGCCTTGGGGAAATGA
AAAGGAAAAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



View online »

Protein Sequence:	>MR212960 representing NM_206823 Red=Cloning site Green=Tags(s) MNIGNCSLNEFIFVGVTNNPEMKGTLFTIFLLIYLINLLGNIGMIILIRMDPQLHTPMYFFLSHLSFCDL CYSTAIGPKMLVDMFGKNKSIPFWGCAVQFFISCTFADSECVLLAVMAFDYQAI SNPLLYTANMSSRKC FMFMAGVYLVGTSDALIHTTLAFRLCFCGSNEINHFFCDLPPIYLLSCSDTQVNYLALFTIYGFLSTI SGVLVSYCYIISSVLKIRSTEGRFAFSTCTSHLTAVAI FQGTLLFTYFRPSSSYSLDQDKMTSLFYTLV IPMLNPLIYSLRNKDVKEALGKMKRKR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_206823
ORF Size:	921 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_206823.1, NP_996559.1</u>
RefSeq Size:	924 bp
RefSeq ORF:	924 bp
Locus ID:	110511
Cytogenetics:	2 D
MW:	34.7 kDa