

## Product datasheet for **MR212973**

### **Olfr944 (NM\_146507) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Olfr944 (NM\_146507) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Olfr944  
**Synonyms:** MOR171-19  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR212973 representing NM\_146507  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGCAGCCGAAACCACTGCACAGTGACTGAGTTCTTCTTAGTTGGGCTCTCAGAGAAGCCAGAATTCC  
AGCTGCCCTCTCCTCCTCTTCTTGAATCTATCTAATCACTGTGACAGGAAACGTGGGCATGATCAC  
ACTGATTGGGCTTAGTCCCACCTGCACACCCCATGACTTTTTCGTCAGGAGTCTGTCCTTCATTGAC  
TTCTGTCAGTCCACAGTTGTTATTCCTAAAATGCTCATGAGCTTCTGACAGAGAAGAACATCATTTCT  
ACTCTGGATGCATGGTTCAGCTCTACTTCTTTCATATTTGGTATTGCAGAGTGCTACAGTTAGCTGC  
AATGGCCTATGACCGATATGTTGCTATTTGTAACCCCTTGCTTTATAATGTAACCATGTCCTATCAGATT  
TACAATCTCTGATTTCCGGGTCATATATTTTGTCTGTGGTCTGTTTCATCCTTAATCACTGGCTTCATG  
TTAGGATTCAGTTCTGTAATTTAGATGTGATTAACCACTATTTCTGTGATCTTCTCCCCTCTTGAATCT  
TGCATCCTCTAATATCTATCAATGAAATATTGATTCTAGTTATTGCTACACTGAATGTCTTTATCCCA  
GTGATGACCATTACTTCTACATCTTATTGCCCACCTCTCTACATTCAGTCCAGTGAGGGCA  
AGTTCAAAGGGTTTAGTACTTGTAGTACCCACATCTCTGCTGTTGCTATCTTTATGGTTTCAGGAGCATT  
CACATACTTACAGCCCTCATTACTGAATTCTATGGGCAAGCAAAAGTGTCTCTGTGTTTTGTACTACA  
GTTGTACCATGCTGAACCCCTTCACTACAGCCTGAGGAATAAGGATGTCAGTTTTGACTGAAAAAAA  
TATTTGAAAGAAAAAATTCATG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR212973 representing NM\_146507  
Red=Cloning site Green=Tags(s)

MAAGNHCTVTEFFLVGLSEKPEFQLPLFLLFLGIYLITVTGNVGMITLIGLSSHLHTPMYFFVRSLSFID  
 FCQSTVVIPKMLMSFLTEKNIISYSGCMVQLYFFFIFGIAECYLAAMAYDRYVAICNPLLYNVMTMSYQI  
 YNSLISGSYIFAVVCSLLITGFMFRIQFCNLDVINHYFCDLLPLLNLASSNIYINEILILVIATLNVFIP  
 VMTIITSYIFIIATILYIHSSEKFKGFCSTSTHISAVAI FYGSGAFTYLP SLLNSMGQAKVSSVFCTT  
 VVPM LNPF IYSLRNKDV SFALKKIFERKKFM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9067\\_f04.zip](https://cdn.origene.com/chromatograms/mm9067_f04.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_146507

**ORF Size:** 933 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_146507.1](#), [NP\\_666718.1](#)

**RefSeq Size:** 936 bp

**RefSeq ORF:** 936 bp


**Locus ID:** 258500

**Cytogenetics:** 9 A5.1

**MW:** 35 kDa

**Gene Summary:** Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

**Product images:**



Circular map for MR212973