

Product datasheet for **MR212978**

Olfr959 (NM_146508) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr959 (NM_146508) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr959
Synonyms:	M34; MOR224-3; ORL583
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR212978 representing NM_146508 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGGAATCGTTCGCTAGTGACTCAGTTTATCCTGCTCGGCATCCCAAACACAGAGGGTCTGGAACTA
 TGCTCTTTGCTGTTTTTGTCTTTCTACATCTTCACTCTGATGGGAACTTATTAATCTTACTGGCAAT
 TATCTCTCCAGTCGGCTTCATACCCCATGTACTTCTCCTATGTAATTGTCTATTTTGTGACATATTT
 TTCCCTTCTGTGAGTTCTCCCAAGATGCTGTTCTACCTCTCAGGAAACAGCAGAGCCATCTCCTATGCAG
 GCTGTGTGAGCCAGCTCTTCTTACCATTTCCTTGGTTGTACCGAGTGTTCCTGTACACAGTGATGGC
 ATATGACCGTTTTGTGGCATATGCTACCCATTACGCTACTCAATAATTATGAGTCACAGAGTTTGTGCG
 ATCCTGGCCACGGGGACATCCTTTTTTGCTGCATTACAGGCTACTTTCTGACCACTCTCACCTTCCAGT
 TGCCCTACTGTGGTCCCAATGAGGTGGACTACTTCTGTGATATCCCTGTGATGCTAAACTGGCTTG
 TGCAGACACATCAGCCCTGGAGATGGTGGGTTTCATCAGTGTGGGCTGATGCCGCTCAGTTGTTTCCTC
 CTCATCTCACCTCCTACAGCTGCATCGTCTGCTCTATTCTGCAGATCCGATCTGCTGAGGGCCGTCGTA
 GAGCCTTCTCCACCTGCAGTGCCACCTCACTGCCATCTTGCTTTTCTACATGCCAGTTGCCTCATATA
 CTTAAGACCAACTCCAAGCCCTTGGCTGGATGCAACTGTGCAGGTCTTAATAATCTGGTCACTCCTATG
 CTAACCCATTAATCTACAGTCTCAGGAATAAAGAGGTGAAAGCATCACTGTGGAAGGTCTACGTAAC
 CTGCCTTTGCTCCTGAGCAGTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR212978 representing NM_146508
 Red=Cloning site Green=Tags(s)


MRNRSVVTQFILLGIPNTEGLETMLFVLFLSFYIFTLMGNLLILLAISSSR LHTPMYFFLCKLSIFDIF
 FPSVSSPKMLFYLSGNSRAISYAGCVSQLFFYHFLGCTECFLYTVMAYDRFVAICYPLRYSIIMSHRVCA
 ILATGTSFFGCIQATFLTTLTFQLPYCGPNEVDYFCDIPVMLKLACADTSALEMVGFISVGLMPLSCFL
 LILTSYSCIVCSILQIRSAEGRRRAFSTCSAHLTAILLFYMPVVLIIYLRPTSPWLDATVQVLNNLVTPM
 LNPLIYSLRNKEVKASLWKVLRKPAFAPEQL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9067_f05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_146508

ORF Size: 933 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_146508.1](#), [NP_666719.1](#)

RefSeq Size: 936 bp

RefSeq ORF: 936 bp


Locus ID: 258501

Cytogenetics: 9 A5.1

MW: 35.2 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR212978