

## Product datasheet for **MR213010**

### Fance (NM\_001163820) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Fance (NM_001163820) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Fance
Synonyms:	2810451D06Rik; AI415634; AW209126
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MR213010 representing NM\_001163820  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

CTGGCGACGCCGGGTTCGGACTCCGTAGCCTCGGCCGAGCGGCCCTGGGCGAGCCTGGAGGCCCTG  
 CCCGCCTCTGCTTCAGGCGCTACAGGCAGGGCCTGAGGGTGCGCGGCGGGGCTAGGGTTCTGCGCGC  
 CCTCGGCCGCGCGGGAGCACTTCCCTTGGGATGGCTTCTGGAGGCGCTGGGCCACTTAGAGCCCGAG  
 GTGCGGGGCCCGACGGCCGCTAGAGCTGGTCCCACTGTTGCTCCGGCTGCCTGGGGTCTGCCAGAAGA  
 ACCTGTGTCCCTGCTGCTGGCTCTCTACCATCTTACCTGAAAGCGGACTCCGCTCTGTGCTGCAGCT  
 CCACCATCAGGATGTGCTCTACCACTGATGCCTGGCTCCATGCCCTGGGGGAGCTGCTGCGAAGGGAT  
 GTGGGGTTGGAGTCGCTGAGGGATCTTCTCCATTGACCAGAAGCTGTACAGTTACAGCTCCGGGACCTGT  
 GTGGGCGGCTGGCCAAGGGGGAGGGGCTGAAACTGGCCCTGGCTCCAGATCTGAACAAGAGGACAG  
 ACTCTCACAGCTTTCGCGGAAACGGACGAAAGAGCCAGAAGAGGCTGCCAGCCCTGAGTCAGAGAGATCC  
 CCTAAAAGGTTCCGGGGCTGTGAGGAGGCGGTGGAAGGAAAGGAACCGGAGGAGAGACCCACGCTGGAGT  
 CGCTGGGATCCCCACCAGATGCAGGAGGCGTGTTCCTGACACTGACGCCAGGCTCCGGAGACTGGCCC  
 TGGCGTGGAGGGTCCCAAGGGTCCCGCTGAGAGTGTGGAGTTGCCCAAAGTTGTCCAGGACCAGGTGCC  
 AGGCTGCAGCTGCTGCTGAAGGCCTTCCAGGAGGGGCTGGAGGGTCCAGGAGAAGCCCTAGTGGACCTGC  
 AGTTTCTTCATGAATGTAGTCCCAGCGAGATGGAGTTGCTATGCAGCGAACTACAGCTACCCAGCTCCC  
 TGACGGAGGTCTCTGCAGCTCTGCAGCCACCTGATGGGTCTCACACCAGCCCTCAGCCTCAGCAATGCC  
 TCTGTGCTGGCCAGGAGCCTCTTCTTGACCGGATCCGCTCCCTGCCGTCTTCTGCCTCCAGGCTTCTCA  
 GAGTGGCCCTCGTCTCCTTCTGTGTAAGTACACCTACGCCATCTGCAGGGCTGTCCCTGTCCCTGGCT  
 CCAGGACCCTCGTGTAGTCTCTGCGCAGACCGAGTTACTGTGTTCCCTCATAAAGGATGATCCCTGGAG  
 TCAGACATGCAGGTCCAGATTTTGGGCGAGTCTGGAGCTGGCCTGGCGAGAAGAGACGTTCTGGTGT  
 TGCAGACGCTCTGGAACGGCAGATCACAGAGCAGCAGAGCCTGGACCTGGCTGTGGCCCTAGAGCCCAA  
 CGCCACCTTCTGAAGAAGGCCCTGCAAGCAGCGCTGAGACATGTGACCCAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**

>MR213010 representing NM\_001163820  
 Red=Cloning site Green=Tags(s)

LATPGSDSVASAGAAPWASLEAPARLLLQALQAGPEGARRGLGVLRALGRRAEHFPWDGFLEALGHLEPE  
 VRGPDGRLELVPLLLRLPGVCQKNLLSLLLALLPSLPESGLRSVLQLHHQDVSSTTDAWLHALGELLRRD  
 VGVGVAEGSSPLTRSCQLQLRDLCGRLQGGRGLKLALAPDPEQEDRLSQLCGKRTKEPEEAAPESERS  
 PKRFRGCEEAVEGKEPEERPTLESLSPPDAGGVLPTDAQAPETGPGVEGPKGPAESVELPKVVQDQVP  
 RLQLLLKAFQEGLEGQEKPLVDLQFLHECSPSEMELLCSELQLPQLPDGGLLQLCSHLMGLTPALSLSNA  
 SVLARSLFLDRIRSLPSSASRLLRVALVSFCVKYTYAICRAVLCPLLQDPRVGPATTELLCSLIKDESLE  
 SDMQVQILGQVLELAWREETFLVLQTLLEQRITEQQSLDLAVALPNATFLKKALQAALRHVTH

**TR**TRPLEQ**KL**ISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**


**ACCN:** NM\_001163820

**ORF Size:** 1452 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001163820.1](#), [NP\\_001157292.1](#)

**RefSeq Size:** 1890 bp

**RefSeq ORF:** 1455 bp


**Locus ID:** 72775

**Cytogenetics:** 17 A3.3

**MW:** 52.7 kDa

**Gene Summary:** This gene encodes the complementation group E subunit of the multimeric Fanconi anemia (FA) nuclear complex composed of proteins encoded by over ten Fanconi anemia complementation (FANCD1) nuclear complex composed of proteins encoded by over ten Fanconi anemia complementation (FANCD1) group genes: FANCA, FANCB, FANCC, FANCD1 (also called BRCA2), FANCD2, FANCE, FANCF, FANCG, FANCI, FANCL, FANCM and FANCN (also called PALB2). The FA complex is necessary for protection against DNA damage. This gene product is required for the nuclear accumulation of FANCC and provides a critical bridge between the FA complex and FANCD2. Defects in the related human gene are a cause of Fanconi anemia, a genetically heterogeneous recessive disorder characterized by cytogenetic instability, hypersensitivity to DNA crosslinking agents, increased chromosomal breakage, and defective DNA repair. Translation of this protein is initiated at a non-AUG (CUG) start codon, which is inferred from the related human gene and the notion that this protein is functionally indispensable. Multiple transcript variants encoding different isoforms have been identified. [provided by RefSeq, Aug 2009]

**Product images:**



Circular map for MR213010