

Product datasheet for **MR213066A1V**

Mouse Lurap1 (NM_026547) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Lurap1 (NM_026547) AAV Particle
Tag: Myc-DDK
Symbol: Lurap1
Synonyms: 1520402A15Rik; Lrp35a
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR213066 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGGGGACCGCGGAGTCCCAGACGCCGGACTTGCGAGACGTGGAGGGCAAGGTGGGCAGGAAGATCC
 CTGAAGGACTGCTCCGCGGGCTTCGAGGCGAGTGCGAGCTGGGAACCTCTGGGGACGTGCTGCTCCCGGG
 GGCGCCACGACCGGCCACGGCCTAGGGGACAAGATCATGGCGCTGAGGATGGAGCTGGCTTACCTGCGA
 GCTATCGATGTGAAGATCCTGCAGCAGCTGGTGACCTTGAACGAGGGTATTGAGGCAGTGCCTGGCTGC
 TGGAGGAGCGGGGTACACTGACGAGCCATTGCAGCAGTCTCACCAGCAGTCAGTACAGTCTGACAGGGGG
 AAGCCCAGGACGCTCAAGGCGGGCAGTTGGGACAGCCTGCCAGACACCAGTTCAACTGACCGGCTGGAC
 AGTGTCTCTATTGGCAGCTTCTTGACACTGTGACCCCCAGAGAGCTGGATGAACAAGGACTCCCTGGGC
 CTTCTGTCTGAGATTGACTGGGCAAAGTTATCCCCACTGAGGACAGGGCCAGGACTGAGGTGGACAT
 GACATCCACTAAGCTAGGGAGCTTAACAGCCACATGGAAGCTCCCGGGGATGGGCTCCAGTGTGGACCA
 CCTGAGCCATCAGAGGATGGCAATGCTAACCAGGGCTTTGAAGCTCATTGGTACTGGGGACAGTGCCAGG
 ATGATGTGACCTTCTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>MR213066 protein sequence Red=Cloning site Green=Tags(s) MEGTAESQTPDLRDVEGKVGRKIQEGLLRGLRGECELGTSGDVLLPGAPSTGHGLGDKIMALRMELAYLR AIDVKILQQLVTLNEGIEAVRWLLEERGTLTSHCSSLTSSQYSLTGGSPGRSRRGSWDSLPTSSDRLD SVSIGSFLLDTVTPRELDEQGLPGPSCPEIDWAKVIPTEDRARTEVDMTSTKLGSALTATWKLPGDGLQCGP PEPSEDGNANQGFEAHWYWGQCQDDVTFL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_026547
ORF Size:	717 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_026547.1, NP_080823.1</u>
RefSeq Size:	1427 bp
RefSeq ORF:	720 bp
Locus ID:	68075
UniProt ID:	<u>Q9D6I9</u>
Cytogenetics:	4 D1
MW:	25.8 kDa