

Product datasheet for **MR213112**

Srp72 (NM_025691) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Srp72 (NM_025691) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Srp72
Synonyms:	72kDa; 5730576P14Rik; A132477; BC019196; C77589
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>MR213112 representing NM_025691
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGC**C

ATGGCGAGCGCGGCAGTGGTGGGTGTCACTCCCTGCGCTGTGGAGTGAGGTGAACCGTTATGGGCAGA
 ACAGCGACTTCACGCGTGCCCTCAAGACTGTCAACAAGATATTGCAGATCAACAAGGATGATGCTACTGC
 CCTGCACTGTAAAGTAGTTTGCCTTATCCAGAATGGAAGTTTCAAAGAAGCCTTGAATGTCATCAATACT
 CACACCAAAGTGTAGCCAATAACTCCCTGTCTTTGAGAAGGCATATTGCGAGTACAGGTTGAACAGAA
 TTGAGAATGCCTTGAAGACAATAGAGAGTGCCACCCAGCAGACAGACAACTGAAGGAGCTTTATGGACA
 AGTGCTGTACCGCTGGAACGCTACGATGAGTGCTTGGCTGTGTACAGAGATCTTGTCCGGAACCCAG
 GACGACTATGATGAGGAGAGGAAAAACAACCTGTACGCGTGTGTCGCGCTCAGAGCAACTGGGAAAAAG
 TGGTTCTGAGAACTTGGGTCTCCAAGAAGGCACACAGAGCTCTGTTACAACGCTGCATGTGCACTGAT
 AGGGCAAGGCCAGCTGACCCAGGCCATGAAAATCCTGCAAAAAGCTGAAGATCTTTGTCGCCGTTTCATT
 TCAGAAGATTCAGATGGGGCTGAGGAGGACCCCAAGGCGAACTGGCCATCATCCATGGTCAGATGGCGT
 ACATCATGCAGCTGCAGGCTCGCACAGAGGAGGCGCTGCAGCTCTACAATCAGATCATCAAGCTGAAGCC
 AACAGATGTGGCACTTCTCGCTGTAATCGCAAATAACATCATTACCATTAAATAAGGACCAAAATGTCTTT
 GACTCCAAGAAAAAGTGAATTAACCAATGCAGAAGGAGTGGAGTTCAAGCTCTCAAAGAGACAGCTGC
 AAGCCATAGAGTTTAACAAAGCTTTGCTTGCCATGTATACAAACCAGGCAGAACAGTGCCGAAAAATAGC
 TGCCAGCTTGCACTCGCAGAGTCCTGAGTATCTCTGCTGTGCTAATCCAAGCCGCCAGCTCTGTCGT
 GAAAAGCAGCATACAAAAGCAATAGAGCTGCTTCAGGAATCTCTGATCAACATCCAGAAAATGCAGCTG
 AAATCAAGCTTACCATGGCACAATTGAAAATTTCCAAGGTAATTTTCCAAGCGTGTCTAATCCTGAG
 AAGTATAGAGGAGCTGAGGCACAAGCCAGGCATGGTATCTGCATTGGTGACCATGTACAGCCACGAGGAA
 GATATTGATAGTGCCATCGAAGTCTTCACACAGGCTATCCAGTGGTACCAAAGCCACCAGCCCAATCCC
 CTGCTCATTTGCTCTTGATAAGAGAAGCTGCAAACTTCAAGCTCAAAATATGGCCGGAAGAAGGAGGCAGT
 CAGTGACCTAGAACAGCTGTGGAACAGAATTCAAAAGATATTACACCCTAGCCCAACTTATTTCTGCT
 TACTCGCTTGTTGATCCAGAGAAAGCAAGGCTCTTAGCAACACCTGCCCTCATCGGATAGCATGTCTC
 TAAAAGTAGATGTTGAGGCGCTTGAATTTCTCTGGTGCTACCTACATTGGAAGAAGGGTGGAAAAGT
 TACTGGAGATAACCAACCAAGGAACAAGGACAGGAGATTTGAAAAGAAGAAGAAGAAAAGAAAGGA
 AAAGTGCCTAAGAATTATGACCCAAAAGTTACCCAGATCCAGAAAGATGGCTACCAATGCGAGAGCGCT
 CTTATTACCGAGGAAGAAAGAAGGGTAAAAAGAAGGACCAGATTGAAAAGGAACACAGGGAGCGACTGC
 AGGAGCCTCCTCTGAACTGGATGCCAGTAAAGCCGTGAGCAGCCACCCACCTCCCGAGACCCGGCAGT
 GCGGCAACAATATCGTCCTCTGCAAGTAACATTGTACCCCAAGACACCAGAAACCTGCAGGCGCCCCAG
 CAACCAAAAAGAAACAGCAACAGAAGAAGAAGAAAGGAGGGAAAAGTGGCTGG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR213112 representing NM_025691
 Red=Cloning site Green=Tags(s)


MASGGSGGVSPALWSEVNRYGQNSDFTRALKTVNKILQINKDDVTALHCKVVCLIQNGSFKEALNVINT
 HTKVLANNLSFEKAYCEYRLNRIENALKTIESATQQTDLKEL YGQVL YRLERYDECLAVYRDL VRNSQ
 DDYDEERKTNLSAVVAAQSNWEKVVPENLGLQEGTHELCYNAACALIGQGQLTQAMKILQKAEDLCRRSF
 SEDSDGAEEDPQAE LAI IHGQMAYIMQLQGRTEEALQLYNQIIKLKPTDVALLAVIANNIITINKDQNVF
 DSKKKVKLTNAEGVEFKLSKRQLQAIEFNKALLAMYTNAEQCRKIAASLQSQSPEYLLPVL IQAAQLCR
 EKQHTKAIELLQEFSDQHPENAAEIKLTMAQLKISQGNISKACILRSIEELRHKPGMVSA LV TMY SHEE
 DIDS AIEVFTQAIQWYQSHQPKSPAHL SLIREAANFKLYGRKKEAVSDLEQLWKQNSKDIHTLAQLISA
 YSLVDPEKAKALSKHLPSSDSMSLKVDVEALENSPGATYIRKKGKVTGDNQPKQGGDLKKKKKKKKG
 KLPKNYDPKVTPDPERWLP MRERSYYRGRKKGKKKDQIGKGTQGATAGASSELDASKAVSSPPTSPRPGS
 AATISSASNIVPPRHQKPAGAPATKKKQKQKKKKGGKSGW

SGP TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9047_g05.zip

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_025691

ORF Size: 2013 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_025691.1, NP_079967.1

RefSeq Size: 3605 bp


RefSeq ORF: 2016 bp

Locus ID: 66661

Cytogenetics: 5 C3.3

MW: 75.1 kDa

Product images:



Circular map for MR213112