

## Product datasheet for **MR213115L3V**

### Hsf2bp (NM\_028902) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hsf2bp (NM_028902) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hsf2bp
Synonyms:	4932437G14Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_028902
ORF Size:	1014 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR213115).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_028902.1</a> , <a href="#">NP_083178.1</a>
RefSeq Size:	2867 bp
RefSeq ORF:	1017 bp
Locus ID:	74377
UniProt ID:	<a href="#">Q9D4G2</a>
Cytogenetics:	17 B1
Gene Summary:	Inhibits BNC1 transcriptional activity during spermatogenesis, probably by sequestering it in the cytoplasm (PubMed:23707421). May be involved in modulating HSF2 activation in testis (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)