

Product datasheet for MR213147

Olfr156 (NM_019474) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr156 (NM_019474) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr156
Synonyms:	mOR37b; MOR262-6; Olfr37b; OR37B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR213147 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAGGAGCCAATCAGTCTACGGTAGCTGAATTTGTCTTGCTGGGGTTGTCAGATCACCCAAAGCTGG
AGAAAACCTTCTTCGTGCTCATCCTGCTGATGTACCTGGTGATCCTGCTGGGCAACGGCTCCTCATCCT
GGTGAGCATCCTCGACTCCACCTGCACACGCCCATGACTTCTCCTGGGGAACCTCCTCTCCTGGAC
ATCTGCTATACTACCTCCTCGATTCCCTTAGTCCTGGATGGTTTTCTCACCCAGGAAGACCATCTCCT
TCTCGGGCTGTGCCGTGCAGATGTTTCTCCTTTGCCATGGGAGCCACGGAGTGTGTCTCCTGGGCAT
GATGGCGTTTGATCGTTATGTGGCCATCTGCAACCCTCTTAGATACCCTGTGGTCATGAACAAGTCTGCC
TATGTGCCATGGCCGTACGCTCCTGGGTGGTGGCGGTGCTAACTCCTTGGTGCAGATCTCCCTTGGCG
TGCAGTTGCCCTTCTGTGGGGACAATGTCATCAATCACTTACCTGTGAGATCCTAGCAGTCTTAAAGCT
AGCCTGTGCTGACATCTCCATCAATGTGATCAGCATGGGGGTGGCCAATGTGATTTTCTGGGGTTCCA
GTTCTGTTTCATCTTTGTCTCCTACATCTTCACTCTCCACCATCCTGAGGATCCCCTCTGCCGAGGGGA
GGAAGAAGGCCTTCTCCACCTGTCTGCCACCTCACCGTGGTACTTGTCTTCTATGGGACCATCCTGTT
CATGTACGGGAAGCCCAAGTCCAAGGACCCATTGGGGCAGACAAACAGGATGTTTCAGACAAGCTCATC
TCCTTATTTTATGGGGTCTGACCCCATGCTGAACCCATCATCTACAGCCTGAGGAACAAGGATGTGA
AGGCCGCTGTGAGGAACCTGGTAGGTCAGAAATGCCTGATTGAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR213147 protein sequence
Red=Cloning site Green=Tags(s)

MEGANQSTVAEFVLLGLSDHPKLEKTFVLLILLMYLVILLGNGVLILVSIILDSHLHTPMYFFLGNLSFLD
 ICYTTSSIPLVLDGFLTPRKTIISFSGCAVQMFLSFAMGATECVLLGMMAFDRYVAICNPLRYPVVMNKSA
 YVPMVAVSSWVAGGANSLVQISLAVQLPFCGDNVINHFTCEILAVLKLACADISINVISMGVANVIFLGV
 VLFIFVSYIFILSTILRIPSAEGRKKAFTCSAHLTVVLVVFYGTILFMYGKPKSKDPLGADKQDVSDKLI
 SLFYGVLPMLNPIIYSLRNKDVKA AVRNLVGQKCLIQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_019474

ORF Size: 957 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_019474.2](#), [NP_062347.2](#)

RefSeq Size: 1099 bp

RefSeq ORF: 957 bp


Locus ID: 29846

Cytogenetics: 4 B1

MW: 34.8 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR213147