

Product datasheet for MR213211A1V

Mouse Gm773 (NM_001033423) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Gm773 (NM_001033423) AAV Particle
Tag: Myc-DDK
Symbol: Gm773
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR213211 representing NM_001033423
 Red=Cloning site Blue=ORF Green=Tags(s)

 TTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGACAATAGGAAGCACTTCTCAGTCTATAACCTTGATGATGGTGAAGAAGAGGATGATGTTCTGGAGT
 ATTCAGCACACTCTACACTTACAGCAAGAAGCCTAGTATTTAGTAATACTGAGAAGGCATTGGTTCCTAC
 AAGAATGGGTGGTACAGATAACAAGAAATGAAGTAAAGAATAGGATGGAACATTGGAGCTCACCTTAGC
 AGTGCTGTTGTTGAAAAGAAAAGGAGGATGGAGTCATATATTAAGATGTTGTCAAAGATTGTGACAAGC
 ATATCAACAACATTGGAAAGTTCATCAAGATGAAGTACGGGAGTTTAAATGAATATACTGAGCAGAT
 TATAACTCTGTTTCAGCAGTGGGATTTTGATATGAAGAATGTTGAAGACCTTGAAGATAACCTAACTAGT
 GCCTTTTGTGAGCAACAGAAGACTTTTCCACAACACATCACGGCAGGTCCACAGCCTGAAAGCACTTA
 AGCAGGCAATGACACGCTCTTAAAGGGTATGGAGAACTTGGACAAAGACAATAGCAGTGTACTTTCTAG
 TATTCACAGTGAGTTTAAAAGCAAATGGCTTTCCTGGAAGAAAATTGACAGATGGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR213211 representing NM_001033423
 Red=Cloning site Green=Tags(s)

 MDNRKHFSVYNLDDGEEEDVLEYSAHSTLTARSLVFSNTEKALVPTRMGTDTRNEVKNRNMENIGAHLS
 SAVVEKKRRMESYIKDVVKDCDKHIKQHWKVHQDEVREFKIEYTEQIITLFQQWDFDMKNVEDLEDNLTS
 AFCQQKTFPQTTSRQVHSLKALKQANDTLKGMENLDKDNSSVLSSIHSEFKKQMAFLERKLTG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2


[View online »](#)

ACCN:	NM_001033423
ORF Size:	618 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001033423.2</u> , <u>NP_001028595.1</u>
RefSeq Size:	951 bp
RefSeq ORF:	621 bp
Locus ID:	331416
Cytogenetics:	X A5
MW:	24 kDa