

## Product datasheet for MR213268

### Krtap4-1 (NM\_001048196) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Krtap4-1 (NM\_001048196) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Krtap4-1  
**Synonyms:** A030010K20Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR213268 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGATCGCC

ATGGTCAACTCTTGTGTGGCTCTGTCTGCTCTGAGGAGAGCTGTGGCCAAGGCTGCTGCCAGCCAGCT  
 GCTGCCAGACCCTGCTGTAGGACCACCTGCTGCCGCCCCAGCTGTGTGTGCCAGTCTGTCAGGCC  
 CAGCTGCTGTGTCTAGCTGCTGCAGACCCAGTTGCTGCAGGCCCCAGTGTGCCAGTCTGTGTGCTGC  
 CAGCCACCTGCTGCCGCCCCAGTTGCTGCATTTCCAGCTGCTGTCGACCTTCTGTGGTAGCTCCAGT  
 GCTGTGGTTCTAGCTGCTGCCGCCCTGCTGCCGCCCTGCTGCAGCCCTGCTGCAGCCCTGCTGCC  
 CCCTGCTGCCGTCCCTGCTGCCGCCCTGTTGCCGCCCTGTTGCTGTCTGAGACCAGTCTGTGGTCAG  
 GTCTGCTGCCAAACCACTTGCTACCGCCCCACCTGTGTCATCTCCACCTGCCCCGCCCATGTGCTGTG  
 CCACACCTGCTGCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR213268 protein sequence  
 Red=Cloning site Green=Tags(s)

MVNSCCGSVCSEESCGQGCCQPSCCQTTCRRTTCRPSCCVSSCCRPSCCVSSCCRPSCCRPQCCQSVCC  
 QPTCCRPSCCISCCRPSCGSSCCGSSCCRPCCRPCCSPCCSPCCRPCCRPCCRPCCCLRPVCGQ  
 VCCQTTYRPTCVISTCPRPMCCATPCCC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



**Cloning Scheme:**


**ACCN:** NM\_001048196

**ORF Size:** 510 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001048196.1](#), [NP\\_001041661.1](#)

**RefSeq Size:** 1011 bp


**RefSeq ORF:** 510 bp

**Locus ID:** 665891

**Cytogenetics:** 11 D

**MW:** 17.8 kDa

Product images:



Circular map for MR213268