

Product datasheet for **MR213317A1V**

Mouse Olfr161 (NM_146860) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr161 (NM_146860) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr161
Synonyms:	MOR131-1; OR912-47M6
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR213317 representing NM_146860 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGAGGCACCAACCAAGTCCAGTGTCTCTGAGTTCCTCTCTGGGACTCTCCAGGCAGCCCCAGCAGC
AACAGCAGCTCATCTTTCTACTCTTCTCATCATGTACCTGGCCACTGTCTTGGGAAACCTGCTCATCAT
CCTGGCCATCAGCACAGACTCCCGCTGCACACCCCATGTACTTCTTCTCAGCAACCTGTCTTTGTG
GATGTATGCTTCTCTCCACCACTGTCCCTAAGGTACTGGCCATTACATACTAAGAAATCAAGCCATTT
CGTTCTCTGGGTGCCTCACACAGCTGTATTTCTCTGTGTGTTTGCTGACATGGACAATTTCTGTCTGGC
TGTGATGGCCTATGACCGATTTGTGGCCATATGCCACCTTTACACTACACAACAAAGATGACCCATCAG
CTTTGTGCCTTTCTTGTGTTGGGTCTGGATGGTAGCCAGTCTGAATGCTCTGTTGCACACACTGCTCG
TGGCTCAACTCTACTTCTGTGGGACAATGTGATCCCCCACTTCTTGTGAAGTGACTCCCCTGCTGAA
ACTCTCTTGTCTCAGACACACATCTCAATGAGTTGATGATTCTTGCTGTTGCAGGGCTGATAATGTTAGCT
CCATTTGTTTGCATCCTTGTCTTATATCCTTATTGCTTGTGCCATCCTGAAAATCTCATCCACAGGAA
GATGGAAAGCCTTCTCTACCTGTGGCTCACACTTGGCTGTTGTGTGCCTTTCTATGGCACTATCATATC
CCTGTATTTCAACCCCTCATCTTCTCACTCAGCTGGGAGGGACATGGCAGCTGCCATGATGTACACAGTG
GTGACCCCATGATGAAACCTTTTCATCTACAGCCTGAGGAACAGGGACATGAAAGGGGCTTTAAGGAAAG
TGCTTACCATGAGATTTATATCTACACAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR213317 representing NM_146860 Red=Cloning site Green=Tags(s) MGGTNQSSVSEFLLLGLSRQPQQQQQLIFLLFLIMYLATVLGNLLIILAISTDSRLHTPMYFFLSNLSFV DVCFSSTTVPKVLAIHILRNQAISFSGCLTQLYFLCVFADMDNLLAVMAYDRFVAICHPLHYTTKMT HQ LCAFLVVGSWMVASLNALLHTLLVAQLYFCGDNVIPHFFCEVTPLLKLSGSDTHLNELMILAVAGLIMLA PFVCILLSYILIAAILKISSTGRWKAFTCGSHLAVVCLFYGTIISLYFNPSSSHSAGRDMAAAMMYTV VTPMMKPFIIYSLRNRDMKGALRKVLTMRFISTQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_146860
ORF Size:	939 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_146860.1, NP_667071.1</u>
RefSeq Size:	942 bp
RefSeq ORF:	942 bp
Locus ID:	258859
Cytogenetics:	16 A1
MW:	34.9 kDa