

Product datasheet for **MR213332A1V**

Mouse Olfr770 (NM_146863) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: Olfr770

Synonyms: MOR114-5

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >MR213332 representing NM_146863
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGGAACCATACAGTAACGACATTCATTCTCCTGGGACTAACAGAGGATCCACAAATCAAAGCCTGC
TTTTGATTTTTCTCCTTCTACCTACTTGTGAATATAACTGGTAACCTTGACAATTATCTTACTTACATT
AATAGATCCTCACCTTAACTCCAATGTACTTTTCTTGCAAACTTTTCTTCTTAGAAATCTTGTTT
ACATCAGCTTGATTCTAGATACTTGATAATCTAGCAACTGGTGACAAGACCATTACCTATGGTGCCT
GTGCCAGTCAAGCATTTTTACCGACCTTTTGGAGTAACCTGAATTTTTCTCCTGGCCACCATGTCCTA
TGACCGCTATGTGGCCATCTGCAACCGCTGCATTACACCACCATCATGAGCACTGCCTGCAGAAGACTC
CTTCTTTGCTGTTGGGTGGCTGGGTGATCATCATCTGACCCTTTAGCCTGAGTCAAAATCTGCAAT
TCTGTGATTCTAACATCATTGATAGCTTTCTTTGTGATGTGTCTCCCTTCTGAAGATTTCTGCTCAGA
CACATGGGTCAATTGAGCAGATGGTTATAGGATGTGCCGTGCTGACCTTCATCACAACCCTTTTTTGTGTG
GTTCTTTCTATGTGTACATAATCAAGACCATCTAAAATTCCTTCTGCCAGCAGAGGAAAAAGGCT
TTTCTACCTGCTCTCCACATGATTGTGGTTTCCATCACTTATGGCAGCTGCATCTTCATCTACGTCAA
ACCATCAGCAAAGGATTCAGTGCTATCAACAAGGGGTAATAGTACTAACTACCTCCATTGCTCCCATG
TTGAACCATTCATTTACACCTTGAGGAACAAGCAAGTCAAACAGGCCTTCAACGATTCAATCAAAAAA
TTGCCTTAGAATGCCAAATGGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:	>MR213332 representing NM_146863 Red=Cloning site Green=Tags(s) MRNHTVTTFILLGLTEDPQIQSLLLIFLLLTYYLLNITGNLTIILLTLIDPHLKTPMYFFLQNFSLFLEILF TSACIPRYLYNLATGDKTITYGACASQAFFTDLFGVTEFFLLATMSYDRYVAICKPLHYTTIMSTACRRL LLCCWVAGVIIILPPFSLSQLQFCDSNIIDSFLCDVSPFLKISCSDTWIEQMVIGCAVLTFITTLFCV VLSYVYIIKTILKFPSAQQRKKAFTCSSHMIVVSITYGSCIFIYVKPSAKDSVAINKGVIVLTTSIAPM LNPFIIYTLRNKQVKQAFNDSIKKIALECQNG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_146863
ORF Size:	933 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_146863.1, NP_667074.1</u>
RefSeq Size:	936 bp
RefSeq ORF:	936 bp
Locus ID:	258862
Cytogenetics:	10 D3
MW:	35 kDa