

Product datasheet for **MR213380A1V**

Mouse Snorc (NM_028473) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Snorc (NM_028473) AAV Particle
Tag: Myc-DDK
Symbol: Snorc
Synonyms: Snorc
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR213380 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCATCTTGCTGGCCCTGCGCGTGGCGCTGCTGCTCATCTCTGGGGTCTGGCCCTGCGGTGCTCA
CAGAGGGTCCCCAGGAGCCGATCCACCCTGTGGAACGAGCCTATTGAATTGCCCTCTGGTGAAGGTCC
CCTGGAGAGCACAAAGCCACAACCAGGAATTTGCTGCTCTGGCCCTCCGTTCCCCACCTCTGCTCCAGCA
CCAGAGGACAGCACCCCTCCTGCCAGGGTGGACCAGGATGGAGGATCGCTTGACCAGGCGCCATTGCAG
CCATCGTGATCGCTGCCCTGCTGGCTACCTGCGTGGTGGTGGCACTCGTGGTCTGCGCGCTAAGGAAGTT
TTCCGCTTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR213380 protein sequence
Red=Cloning site Green=Tags(s)

MASCLALRVALLLISGVLAPAVLTEGPQEPDPTLWNEPIELPSGEGPLESTSHNQEFVSGPPFPPTSAPA
PEDSTPPARVDQDGGSLGPGAIAAIVIAALLATCVVLALVVVALRKFSAS

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_028473



[View online »](#)

ORF Size:	360 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_028473.1, NP_082749.1</u>
RefSeq Size:	462 bp
RefSeq ORF:	366 bp
Locus ID:	73234
UniProt ID:	<u>Q9CXL7</u>
Cytogenetics:	1 D
MW:	12.1 kDa