

Product datasheet for **MR213385A1V**

Mouse Ccdc192 (NM_029290) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Ccdc192 (NM_029290) AAV Particle
Tag:	Myc-DDK
Symbol:	Ccdc192
Synonyms:	4930583I06Rik
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR213385 representing NM_029290 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGGTTAGTCGGTGTCTGTCCATGGACAGAGGATCTCAATGGGCTGCCAAGATGGGAGGATGTCACA
GCAAGAAAGTTGTAATCCCAGATATTGAGACCTCAGCGAGATGCAGGTCGATGCTTGGAAAGTTATCAGTC
ATATATTCAGTCACAAAACAGAACATCGGGAACTCACCAGACAGTGGACAGACAGTATTTTCTCTGGAA
CACTTGGAGATTTGCCGAAGGAAGCAGAGGAAAAGGCTAGAGCTCTGTTGGAACAGCTTACAGCTTCTG
AAGCAACAAAATCACAACCTCTTGAGAAGGCTCCATGCTGGAGGGAAGGCTAGAGGATGTGAACCGAAA
GAACGTTGGTGGAGAGTTGTATGAAAACATGGTTCTTGAGAAAGACAAGTGCATTGAGAAGTTGCAAGCT
GAAGTCAAAGCCTCCCAGGAGAACTTAACATCCATAAACTAAAGCACAAAAAGGCAGTAAAAAGCTAC
AAACTGATTTGGCAATAGCTAAGCAAGAGGCTGCCATCACTGTCTGGAAGTCAATGAGAAGATAATGAC
ACTCTGTGAAGGAAGGCCCTTCCCTAGAGCAACCAACTCAGTGGAAAGTACTCTCCGGAGGCTCCCACCT
GTGGAAGAAAGCGACAGGAAGTTTTCCCTGATCATGGAGCTGTCCACCAACTTTCCCTCAGACGGAGA
AGATCACTCAGCTCGAAGACACCTTAGCGGAAAAGGAAAAGAAGATCGAGGAGCTGGAAGCAGAGTGCAC
TTCCAGGTCCTCCCAAGGTGAAGACAGACTTACAGAAGATTTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR213385 representing NM_029290 Red=Cloning site Green=Tags(s) MRLVGVCPTEDLNGPAKMGGCHSKKVVIPDIETSARCRSMLGSYQSYIQSQNRSTSGNSPDSGQTVFSLE HLEICLKEAEEKARALLEQLTASEATKSQLEKVSMLGRLDVDNRKNVGGELYENMVLEKDKCIEKLQA EVKASQEKLNIHKLKHKKAVKKLQTDLAIKQEAAITVLELNEKIMTLCEGRPFPRATNSVEVLSGGLPP VEESDRKFSLIMELSTQLSLQTEKITQLEDTLAEKEKKIEELEAECTSQVPQGEDRLTEDL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_029290
ORF Size:	813 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_029290.1 , NP_083566.1
RefSeq Size:	967 bp
RefSeq ORF:	816 bp
Locus ID:	75444
Cytogenetics:	18 D3
MW:	30.2 kDa