

Product datasheet for MR213439

Olfr1297 (NM_146888) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr1297 (NM_146888) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr1297
Synonyms:	MOR248-4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR213439 representing NM_146888 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGATGGAGACAATCAGACTGTTGTGTCAGAATTTGTACTCTTGGGACTTTCGAACTCTAAGAATCTTC
 AAGTCTTACTCTTCTCATCTTCTATTGCTTTATCTGCTCATCATGTCTGGAAATATTGTCATCCAAAT
 CTTAATAACAACCTGACCCCATCTGCATTCTCCCATGTACTTCTTGTGGCCAACTGTCCTTTGTTGAT
 ATGTTGCTCTCCTCAAACACAACCTCCTAGGATGATCATAGATTTTTTCAGGGAAAAGAAGACCATTTCCT
 TTGCAGGATGCATGTCCAGATCTTTTTTCCCACTGCATTGCTTCAGGAGAGGTGGTATTGTTGGCTTT
 AATGGCCTATGACCGCTATGTGGCCATCTGCAAACTCCATTACTTCACCATCATGAACCTGAAGAGG
 TGCAGTGGGTGGTGTGACTTCTTGGACAATTGGCTTTTACATGGTATAAGTCACGTGGTAGTGCTTT
 TGCAGTTCCTTTCTGTGGTCCCAACAAAATAGACAGTTTCTTTTGTGACATGCCCTGGTAATCAAGCT
 AGCCTGCATGGATTCCCAGGATTTAAATACTTTAATGAATGGTGAATGTGGTATTTTGGCTGTAACCTGT
 TTTATACTGTTGCTTATTTCTACACATATATCCTTATCACTGTTTCATCAGAATTCGAAACTGGAGCTT
 CCAAGGCTCTGTCCACATGCACTGCACACATCACAGTGGTAATGATCTTCTTTTGCCTGTTTCTTCAT
 CTATGTGTTTCCCTAAATATCACCTGGTTGGACAAATTTCTTGTGTGTTTATTCTGTTATTACACCT
 CTCCTAAACCCAGTCATTTATACACTGAGAAATAAGAGATGAAAAATGCTATGAAAAGATTCAATTGGCA
 AATTTCTGAGAGCCAAAGGAAATTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR213439 representing NM_146888
 Red=Cloning site Green=Tags(s)


MDGDNQTVVSEFVLLGLSNSKNLQVLLFLIFLLLYLLIMSGNIVIQILITDPLHSPMYFLLANLSFVD
 MLLSSNTTPRMIDFFREKKTISFAGCMSQIFFSHCIASGEVLLALMAYDRYVAICKPLHYFTIMNLKR
 CTGLVLTSWTIGFLHGISHVVLLQLPFCGPNKIDSFFCDMPLVIKLACMDSQDLNLMNGECGILAVTC
 FILLISYTYILITVHQNSKTGASKALSTCTAHITVVMIFFLPCFFIYVFLNITWLDKFLAVFYSVITP
 LLNPVIYTLRNKEMKNAMKRFIGKFLRAKGNS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9069_d09.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_146888

ORF Size: 936 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_146888.1](#), [NP_667099.1](#)

RefSeq Size: 939 bp

RefSeq ORF: 939 bp


Locus ID: 258890

Cytogenetics: 2 E3

MW: 35.2 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR213439