

Product datasheet for **MR213439A1V**

Mouse Olfr1297 (NM_146888) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr1297 (NM_146888) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr1297
Synonyms:	MOR248-4
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR213439 representing NM_146888 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGATGGAGACAATCAGACTGTTGTGTCAGAAATTTGACTCTTGGGACTTTGAACTCTAAGAATCTTC
AAGTCTTACTCTTTCTCATCTTCCTATTGCTTTATCTGCTCATCATGTCTGAAATATTGTCATCCAAAT
CTTAATAACAACCTGACCCCATCTGCATTCTCCCATGTACTTCTTGTGGCCAACTGTCCTTTGTTGAT
ATGTTGCTCTCCTCAACACAACCTCTAGGATGATCATAGATTTTTTCAGGGAAAAGAAGACCATTTCCT
TTGCAGGATGCATGTCCAGATCTTTTTTCCCACTGCATTGCTTCAGGAGAGGTGGTATTGTTGGCTTT
AATGGCCTATGACCGCTATGTGGCCATCTGCAACCACTCCATTACTTCACCATCATGAACCTGAAGAGG
TGCACCTGGGTTGGTGTGACTTCTTGGACAATTGGCTTTTTACATGGTATAAGTCACGTGGTAGTGCTTT
TGCAGCTTCCTTTCTGTGGTCCCAACAAAATAGACAGTTTCTTTGTGACATGCCCTGGTAATCAAGCT
AGCCTGCATGGATTCCAGGATTTAAATACTTTAATGAATGGTGAATGTGGTATTTTGGCTGTAACCTGT
TTTATACTGTTGCTTATTTCTACACATATATCCTTATCACTGTTTCATCAGAATTCGAAAACCTGGAGCTT
CCAAGGCTCTGTCCACATGCACTGCACACATCACAGTGGTAATGATCTTCTTTTGGCCTGTTTCTTCAT
CTATGTGTTTCCCTAAATATCACCTGGTTGGACAAATTTCTGCTGTGTTTTATTCTGTATTACACCT
CTCCTAAACCCAGTCATTTATACACTGAGAAATAAGAGATGAAAAATGCTATGAAAAGATTATTGGCA
AATTTCTGAGAGCCAAAGGAAATTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR213439 representing NM_146888 Red=Cloning site Green=Tags(s) MDGDNQTVVSEFVLLGLSNSKNLQVLLFLIFLLLYLLIMSGNIVIQILITTDPHLHSPMYFLLANLSFVD MLLSSNTTPRMIDFFREKKTISFAGCMSQIFFSHCIASGEVLLALMAYDRYVAICKPLHYFTIMNLKR CTGLVLTSWTIGFLHGISHVVLLQLPFCGPNKIDSFFCDMPLVIKLACMDSQDLNTLMNGECGILAVTC FILLISYTYILITVHQNSKTGASKALSTCTAHITVVMIFFLPCFFIYVFPLNITWLDKFLAVFYSVITP LLNPVIYTLRNKEMKNAMKRFIGKFLRAKGNS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_146888
ORF Size:	936 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_146888.1, NP_667099.1</u>
RefSeq Size:	939 bp
RefSeq ORF:	939 bp
Locus ID:	258890
Cytogenetics:	2 E3
MW:	35.2 kDa