

Product datasheet for **MR213466A1V**

Mouse Zfp703 (NM_001110508) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Mouse Zfp703 (NM_001110508) AAV Particle

Tag: Myc-DDK

Symbol: Zfp703

Synonyms: 1110032O19Rik; AI430822; AL022941; Csmn1; End2; Zeppo1; Znf703; Zpo1

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >MR213466 representing NM_001110508
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGCGATTGCGCCGCTGGATCTAACCCAAGGACACCCGAAAGCAGCGGCAGCGGCGGCGGCAGCAGCA
GCGGCGGCGGCGGCGGAAGAGGCCGCGGTGCCGGCAGTGGTGTCCCTCCTGCCCGCGGACCCCT
GCGCCAGGCGAACCGCTCCCTATCAGGGTCTGAAGATGCTGAGCGCTCACACCGGTCACCTCCTGCAC
CCGAGTACCTGCAACCGCTGCTCTCCACCCGTCAGCCCAATTGAGATGGGAAATTATCTGGAAAGG
GGCGAATCACCATCCGATTACAGAGAAGTCACCTTTAACCGGAGGGCTGTGGTTGAGCCATTTCTCT
GCGGCCAGCGCCCATGTCACTCGGGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR213466 representing NM_001110508
Red=Cloning site Green=Tags(s)

MSDSPAGSNPRTPESSGSGSSSGGGGKRPAPVAVVSLPPADPLRQANRLPIRVLMKLSAHTGHLLH
PEYLQPLSSTPVSPITEMGKLSGKGRITIRITEKSPLTGRAVVEPFLLRASAHVTLGD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse

Serotype: AAV-2

ACCN: NM_001110508



[View online »](#)

ORF Size:	381 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001110508.1, NP_001103978.1</u>
RefSeq Size:	3320 bp
RefSeq ORF:	384 bp
Locus ID:	353310
Cytogenetics:	8 A2
MW:	13.5 kDa