

Product datasheet for **MR213468**

Olfr884 (NM_001011798) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr884 (NM_001011798) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr884
Synonyms:	MOR162-9P; MOR162-11P; MOR162-13; Olfr1550-ps1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR213468 representing NM_001011798 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTTTAAGAAATGCCTCTTCAGTGAAAGAATTTATCTTGCTGGGATTGACACAGCAGCCAGGGCTCC
 AGCTGCCGCTCTTCTCCTGTTCTTGGGAATCTATGTGGTCTCCATGTTGGGAACCTGGGCTTGATTGT
 TCTGATTGTGTTGAATCCTCACCTGCACACCCCATGTACTACTTTCTTCAACCTTTCCTTCATAGAT
 CTCTGCTACTCCTCTGTCATAACCCCTAGAATGTTGGTGGGTTTTGTGAAGCAGAACATCATCTCTCATG
 CTGAGTGCTTGACTCAGCTTTTTTCTTTGCCTCTTTGTTATTGATGAATGCTACATTTTGACAGCAAT
 GGCTTATGACAGATATGCTGCCATTTGTAAGCCCTGCTTACCAGGTCACCATGTCTCATCAGGTCTGC
 CTATTGATGACTATGGGTGTGTATGTGATGGGCTTTCAGGTGCCTTGTCACACATAGTTTGATGCTGA
 GACTCACCTTCTGTGATGGTAACATCATTAATCACTATGTTTGTGATGTACTTCTCTCTTAACTCTC
 CTGCACAAGTACCTCCATCAATGAGATGGTAGTTTTATTGTTGTGGGTGTCAATGTGATAGTGCCAGC
 CTGACTCTCTTTGTTTCTTATACCTTAATCCTTCCAACATCCTCAGCATCCATTCTGCAGAAGGTAGAT
 CAAAAGCCTTCAGTACCTGTGGCTCCCATGTGATGGCTGTTTCTTTTCTTTGGAGCTGCAGCCTTCAT
 GTATCTTAAGCCTTCTAGTGCATCTGTGGATGAAGAGAAATTATCTACCATCTTTATACCATTTTGGGT
 CCAATGCTGAATCCTTTCATCTACAGTATAAGGAATAAGGATGTCCATCTTGCACTGAGAAAAACATTGA
 TGAAGTGAAGTTTTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR213468 representing NM_001011798
 Red=Cloning site Green=Tags(s)


MALRNASSVKEFILLGLTQQPGLQLPLFFLFLGIYVVSMLGNLGLIVLIVLNPHLHTPMYYFLFNLSFID
 LCYSSVITPRMLVGFVKQNIISHAECLQLFFFAFFVIDECYILTAMAYDRYAAICKPLLYQVTMSHQVC
 LLMTMGVYVMGFAGALSHIVCMRLRTFCDGNIINHVCVDVLPPLKLSCTST SINEMVVFIVGVNVIVPS
 LTLFVSYTLILSNILSIHSAEGRSKAFSTCGSHVMAVSFFFGAAAFMYLKPSSASVDEEKLSTIFYTILG
 PMLNPF IYSIRNKDVHLALRKTLMKLRFS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9066_a10.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001011798

ORF Size: 927 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001011798.1](#), [NP_001011798.1](#)

RefSeq Size: 930 bp

RefSeq ORF: 930 bp


Locus ID: 257996

Cytogenetics: 9 A4

MW: 34.6 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR213468