

Product datasheet for **MR213468A1V**

Mouse Olfr884 (NM_001011798) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr884 (NM_001011798) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr884
Synonyms:	MOR162-9P; MOR162-11P; MOR162-13; Olfr1550-ps1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR213468 representing NM_001011798 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCTTTAAGAAATGCCTCTTCAGTGAAAGAATTTATCTTGCTGGGATTGACACAGCAGCCAGGGCTCC
AGCTGCCGCTCTTCTTCTGTTCTTGGAATCTATGTGGTCTCCATGTTGGGAACCTGGGCTTGATTGT
TCTGATTGTGTTGAATCCTCACCTGCACACCCCATGTACTACTTTCTTCAACCTTTCCTTCATAGAT
CTCTGCTACTCCTCTGTCATAACCCCTAGAATGTTGGTGGGTTTTGTGAAGCAGAACATCATCTCTCATG
CTGAGTGCTTGACTCAGCTTTTTTCTTTGCCCTTTGTTATTGATGAATGCTACATTTTGACAGCAAT
GGCTTATGACAGATATGCTGCCATTTGTAAGCCCTGCTTTACCAGGTCACCATGTCTCATCAGGTCTGC
CTATTGATGACTATGGGTGTGTATGTGATGGGCTTTGCAGGTGCCTTGTCACACATAGTTTGCATGCTGA
GACTCACCTTCTGTGATGGTAACATCATTAACTACTATGTTTGTGATGTAATCTCTCTTAACTCTC
CTGCACAAGTACCTCCATCAATGAGATGGTAGTTTTATTGTTGTGGGTGTCAATGTGATAGTCCCAGC
CTGACTCTCTTTGTTTCTTATACCTTAATCCTTTCCAACATCCTCAGCATCCATTCTGCAGAAGGTAGAT
CAAAAGCCTTCAGTACCTGTGGCTCCCATGTGATGGCTGTTTCTTTTCTTTGGAGCTGCAGCCTTCAT
GTATCTTAAGCCTTCTAGTGCATCTGTGGATGAAGAGAAATTATCTACCATCTTTTATACCATTTTGGGT
CCAATGCTGAATCCTTTCATCTACAGTATAAGGAATAAGGATGTCCATCTTGACTGAGAAAAACATTGA
TGAAGTGAAGTTTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR213468 representing NM_001011798 Red=Cloning site Green=Tags(s) MALRNASSVKEFILLGLTQQPGLQLPLFFLFLGIYVVSMLGNLGLIVLIVLNPHLHTPMYYFLFNLFSFID LCYSSVITPRMLVGFVKQNIISHAECLTQLFFFAFFVIDECYILTAMAYDRYAAICKPLLYQVTMSHQVC LLMTMGVYVMGFAGALSHIVCMLRLTFCDGNIINHYVCDVLP LLKL SCTST SINEMVVFIVVG VNVIVPS LTLFVSYTLILSNILSIHSAEGRSKAFSTCGSHVMAVSFFFGAAAFMYLKPSSASVDEEKLSTIFYTILG PMLNPFYISIRNKDVHLALRKTLMKLRFS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001011798
ORF Size:	927 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001011798.1, NP_001011798.1</u>
RefSeq Size:	930 bp
RefSeq ORF:	930 bp
Locus ID:	257996
Cytogenetics:	9 A4
MW:	34.6 kDa