

Product datasheet for **MR213527**

Frrs1l (NM_001142965) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Frrs1l (NM_001142965) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Frrs1l
Synonyms:	6430704M03Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR213527 representing NM_001142965 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGGGCAACCCCTGCGGCGCCCGGCTGGGTACCACTGCTCCTGCGACTGCTGCTGGCGGGGATCG
CTGCCTGCGATGCCAGCCCCGCGGACGACAGCGCGGGCCCGGGGGCCGGGGACCCCGGGACGCGCGCG
AGGGGACGCGGGCGGACGAGGCAGTGCCGCGCCAGACTCCTCTACGGCACCTTCGCCAGCGAGTTC
TACGACCTGCGCTACCTGTCAGAGGAGGGTTACCCTTTCCCACTGCTCCTCTGTGGACCCATTGCCA
AAATCAAGGTGGAAGACTGTGGAAGAACTAAGGGATGCTTCAGATACGGGAAGCCAGGCTGTAATGCAGA
AACCTGTGACTACTTCCTTAGCTACCGGATGATAGGCGCTGACGTGGAGTTTGAGCTGAGTGCAGACACA
GATGGTTGGGTCGCGGTTGGATTCTCTCAGACAAAAAGATGGGTGGTGACGATGTCATGGCCTGCGTCC
ATGATGACAATGGCAGAGTCCGCATACAGCACTTCTATAATGTGGGTCACTGGGCAAAGGAGTTTCAGAG
AAATCCTGCCCGAGATGAAGAAGGAGTTTTGAGAACAAACCGCTCACCTGCAGGTTTAAAGCGTCCCGTG
AATGTTCCAGGGATGAAACGATCGTGGATCTGCACCTGAGTTGGTACTACCTGTTTGCTGGGGTCCGG
CCATTACAGGGCGCCATCACTCGTCACGACATAGACTCTCCACGGCTTCAGAGCGCGTTGTCAGTATTTA
CAAGTACGAAGACATTTTCATGCCGTGGCTGCCTACCAGACCTCTCCTCGCCATTTGTCTGCTTCTC
ATTGTGCCCTGACCTTCTACCTGCTGATGGGAACCCCT

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR213527 representing NM_001142965
 Red=Cloning site Green=Tags(s)


MAGQPLRRPAWVPLLLRLLLAGIAACDASPADDASAGPGGRGPRGRARGDAGADEAVPRHDSSYGTFASEF
 YDLRYLSEEGYPFPTAPPVDPFAKIKVEDCGRTKGCFRYKPGCNAETCDYFLSYRMIGADVEFELSADT
 DGWVAVGFSSDKKMGDDVMACVHDDNGRVRIQHFIYVVGQWAKEVQRNPARDEEGVFENNRTVCRFKRPV
 NVPRDETIVDLHLSWYYLFAWGPAIQGAITRHDIDSPPASERVVSIYKYEDIFMPAAAYQTFSSPFCLLL
 IVALTFYLLMGTP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9092_c07.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001142965


ORF Size: 879 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001142965.1, NP_001136437.1</u>
RefSeq Size:	3854 bp
RefSeq ORF:	882 bp
Locus ID:	230235
UniProt ID:	<u>B1AXV0</u>
Cytogenetics:	4 B3
MW:	32.5 kDa
Gene Summary:	Important modulator of glutamate signaling pathway.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR213527