

Product datasheet for **MR213542**

Zfp874b (NM_001076791) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zfp874b (NM_001076791) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Zfp874b
Synonyms:	9630025I21Rik; Rslcan16
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR213542 representing NM_001076791
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGGAAATGCTGTCCTTCCGGGATGTGCCATTGATTTTTCTCCAGAAGAGAGGGACTACCTGGGTC
 CTGCTCAGTGGGATCTGTACAGAGACGTGATGCTGGAGAATTATAGCCACCTTGTGTTCTGGGTCTTGC
 TGTGCTAAGCCATACTTGGTAACATTTCTGGAGCAAAACCAAGGTCTTCAGGTGTGAAAAGCCACGCA
 GCAGCCACCATCCAGGAACAACAGGCAATGAATTTAACAATCACAATGAGGCCATTGATCAGAGCTCTC
 TAAGTATGCAATGCCAGAGAATTCATACAGGAGAGGAACCTACTACTTTGAAGACTGTGGTAAGGCCCT
 TAGTTCTCAGGCAACACTTTCTGTACTCCAGAGACGTTATAGTGGAGACAAACCTACAAGTGTAAAGAA
 TGTCAAAAACACTTAGCAGTCGCTCATCACTTTTAATACACCAGAAATATCATACTGATGAGAAAACCT
 ACAAGTGTGAAAAATGTGGCAAACGGTTTTTGCCTTCTCACATCTTCAGCATCATCAAAAAATTCATAC
 TGGAGAGAAACCATACAAATGTGAAGAGTGTGACAAAGCCTTTCTTCATCACTCATATCTTAGGAAACAC
 CAGGCAGTTCATACTGGAGAGAAACCTACAAGTGTGAAGAGTGTGGAATTCATTTTACTATCCTGCAA
 TGCTGAAGCAACATCAAAGAATTCATTCTGGAGAGAACTTGACAAGTGTGAAGAATGTGGGAAAGTATT
 TTCCTCTGCTTTCTTTCTTAACCAGCATAAAGGAATTGATTCCATAGAGAAAAGGTACAAGTGTCAAGAA
 TGTGGCAAATCTTTTGTATCGTTTCATATCTTAGGGAACATTACAGAAGGCATTCTGGAGAGTATCCCT
 ACAAGTGCAGAAGATGTGGCAAAGGTTTTCCCGTTCTCAAAGCTTCAGGAACATCAAACAATTCATAC
 TGGAGTTAAACCTTACAAGTGTGAAGAGTGTGAGAAATGTTTTCTCTTTACGTCCCTTAAAGACAT
 CAAATCATCCATTCTGAAGATACTCCCCATGAGTGTGTAGAATGCGGCAAAAGGTTTTCTAGTTCTTCCC
 GCCTTCAGGAACATCAGAAAATTCACACTGAAGAGAAACCATACAAATGTGAGGAATGTGACAAAGCCTT
 TCTTTATCACTCGTTTCTTAGGAGACACAGAGCAGTTCATACTAGAGAGAAACCTACACGTGTGAAGAG
 TGTGGGAAATGTTATACCTCTTTCACATCCCTGAAAAGACATCAAACAATTCATTCTGTAGACACTCCCC
 ATGAGTGGGTA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR213542 representing NM_001076791
 Red=Cloning site Green=Tags(s)


MEEMLSFRDVAIDFSPEERDYLGAQWDLYRDVMLENYSHLVFLGLAVAKPYLVTFLEQNQGSSGVKSHA
 AATIPGTTGNEFNNHNEAIDHSSLSMQCQRIHTGEEPPYFEDCGKALSSQATLSVLQRRYSGDKPYKCKE
 CHKTLSSRSSLLIHQKYHTDEKTYKCEKCGKRFLRSSHLQHHQKIHTGEKPYKCEECDKAFLHHSYLRKH
 QAVHTGEKPYKCEECGNSFYYPAMLKQHRIHSGEKLKCEECGKVFSSAFFLNQHKGIDSIEKRYKCQE
 CGKSFCYRSYLREHYRRHSGEYPYKCEECGKGFSSSKLQEHQTIHTGVKPYKCEECCKCFSSFTSLKRH
 QIIHSEDTPHECVCGRFSSSSRLQEHQKIHTEKPYKCEECDKAFLYHSFLRRHRAVHTREKPYTCEE
 CGKCYTSFTSLKRHTIHSVDTPEHWV

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI


Cloning Scheme:



ACCN:	NM_001076791
ORF Size:	1341 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001076791.2 , NP_001070259.1
RefSeq Size:	3792 bp
RefSeq ORF:	1344 bp
Locus ID:	408067

Cytogenetics: 13 B3
MW: 52.3 kDa

Product images:



Circular map for MR213542