

Product datasheet for **MR213605L3V**

Mgst2 (NM_174995) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mgst2 (NM_174995) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mgst2
Synonyms:	GST2; MGST-II
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_174995
ORF Size:	441 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR213605).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_174995.3 , NP_778160.2
RefSeq Size:	614 bp
RefSeq ORF:	444 bp
Locus ID:	211666
UniProt ID:	A2RST1
Cytogenetics:	3 C



[View online »](#)

Gene Summary:

Catalyzes several different glutathione-dependent reactions. Catalyzes the glutathione-dependent reduction of lipid hydroperoxides, such as 5-HPETE. Has glutathione transferase activity, toward xenobiotic electrophiles, such as 1-chloro-2, 4-dinitrobenzene (CDNB). Catalyzes also the conjugation of leukotriene A4 with reduced glutathione to form leukotriene C4 (LTC4) (By similarity). Involved in oxidative DNA damage induced by ER stress and anticancer agents by activating LTC4 biosynthetic machinery in nonimmune cells (PubMed:26656251).[UniProtKB/Swiss-Prot Function]