

Product datasheet for MR213635

Olfr1162 (NM_001011835) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Olfr1162 (NM_001011835) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Olfr1162
Synonyms: MOR174-14
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR213635 representing NM_001011835
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTGCCAAATGGAAAGAAATGTGAGTGTGGAGATTATCTTCTCCCTGTTGGGCTTCACAGACTACCCTG
AGCTTCAGATCCCCCTTTTCTAATCTTCTTTCATGTACGTTATCACTGTAGTGGAAACCTGGGAT
GATAGCACTCATCAATATCAACCCTAAGTTCACACTCCAATGTACTTCTCCTTAGCCACCTTTCCTT
GTTGACTTCTGCTACTCAACCATCATCACACCAAGCTGCTTGAAAACCTGGTCTGGCTGATAAACTA
TACTTTACTTCAGCTGCATGTTTCAGTACTTCCCTGTCTGTGGCTGTGGTGTCTGAGTCTACTTACT
AGCAGTTATGGCTTATGATCGCTTCGTGGCCATCTGTAATCCCTTGTGTATACTGTGGCCATGTCACCA
AGGCTGTGCATCCTGCTTGTGACTGGGTATATCTGGAGCATATTTGCTCCATTGATTCTCCTTTGCT
ATGCTCTGCAGCTAAAGTTTTCAAGATTTAATGTGATCAACCACTTCTTCTGTGAATATACTGCCCTCAT
TGCTGTCTCTAGTTCTGATATACACATTTCCAGCCTGCTACTCTTTTGTGTTTGCAACCTTCAATGAAGTG
AGTACACTCTTGATCATTCTCACATCTTATGTGTTCAATTTTCGTGACGGTATTAATAAATCAAATCTGCCA
GTGGACGTCGCAAAGCCTTCTCTACCTGTGCTTCCACCTGACTGCCATCACAATTTTTCACGGAACCAT
ACTTTCTTATACTGTATACCTAACTCCAAAACTCTAGGCAAGTCGTTAAAGTGGCCTCTGTCTTTTAC
ACTGTTGTCAACCCATGTTGAACCCACTGATCTACAGCCTAAGAAACAAGGATGTGAAGGATGCTTTCC
AGAAATTAGTGAGCACAAGATCCCTCTCAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001011835.1](#), [NP_001011835.1](#)

RefSeq Size: 945 bp

RefSeq ORF: 945 bp


Locus ID: 258105

Cytogenetics: 2 D

MW: 36.1 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR213635