

Product datasheet for MR213660

Saysd1 (NM_026209) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Saysd1 (NM_026209) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Saysd1
Synonyms: 1810063B07Rik; 4930488P03Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR213660 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAACAGCGGTTAGCCGAGTTCGGGAGGCCGCAACGGGCCAGTCTGGTAGCCAGCCTTCCACGT
 CCAGTCAGAGCGTACAGACCTCGGGAGCGAAGGCGGAGCCAGCGGCGGCACCCAAAGACAGCCACGGG
 CTGGCTGACACGCTTCTGAAACGGAAAGCCAACCCGCCATCGCTCAGGCTCAGCCAACAGCCTCAG
 GAAGCAGGCCAGCAGCTCCCTGAGAGCACAGCAGTCCCTCTGCCTTCTCCTGCCCCAGTCTTCTCTGA
 CTAACATCACCTTCTTGAAGTTCTCCTCTGGTTGGTCTGCTGGGACTGTTTGTGGAAGTGAATTTGG
 TCTGGCTTATTTGTCTGTCCATGTTCTATTGGATGTATGTAGGACTCGGGGTCCCGAGGAGAAGAAA
 GAAGGGGAGAAGAGCGCCTACTCTGTGTTCAACCCGGGCTGTGAAGCCATCCAGGGCAGCCTGACTGCTG
 AGCAGCTGGAGCAAGAGTTACAGCTCAGACCCCTCAGGGAAGTAGGACCAGCCCACTGTTCTCTCTTA
 TCCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR213660 protein sequence
 Red=Cloning site Green=Tags(s)

MEQRLAEFREARKRASLVAQPSTSSQSVQTSQAKAEPAATPKTATGWLTRFLKRKANPAIAQAQPNQPQ
 EAGQQLPESTAVPLPSSCRQSFLTNITFLKVLWVLLGLFVELEFGLAYFVLSMFYWMYVGTGRGPEEKK
 EGEKSAYSVFNPGEAIIQGLTAEQLELQLRPPQGSRTSPSCSSYP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM_026209

ORF Size: 567 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_026209.1](#), [NP_080485.1](#)

RefSeq Size: 2314 bp


RefSeq ORF: 567 bp

Locus ID: 67509

UniProt ID: [Q8K190](#)

Cytogenetics: 14 A3
 MW: 20.7 kDa

Product images:



Circular map for MR213660