

## Product datasheet for MR213665

### Olfr354 (NM\_146939) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr354 (NM_146939) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr354
Synonyms:	MOR127-4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR213665 representing NM_146939 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGGAAAAATAAGTAGAGTGAACCAATCTGTTGCTTCAGACTTCCTTCTTCTGGGACTCTCTGAACAGC  
CAGGGGAGCAGCCTCTTCTATTTGGCATCTTCTGGGCATGTACCTGGTCACTATGGTGGGGAACCTACT  
CATCATCTTTGCATCAGTTCTGATGCACACCTACATACTCTATGTACTTCTTCTAGCAAATCTCTCC  
TTAACTGATGCCTGTTTCACTTCTGCATCAGTACCCAAAATGCTTGCCAACATTTATACCCAGAGTCAAA  
CCATCTCCTATTCTGGGTGCTTACTCAGTTGTATTTTCTCCTGATGTTTGGCGGCTCGACAACCTGCCT  
TCTGGCTGTGATGGCATATGACCGCTATGTAGCCATCTGCCAGCCACTCCATTACAGCACAGCCATGAGT  
CCCCAACTTTGTGATTAATGCTGTGTACATGCTGGGTGCTGACCAACTGTCTGCCCTTATGCACACAC  
TGCTGCTGACCCGTGGCTTTCTGTGCTCACACAGCCATCCCCATTTCTACTGTGATCCCAGTGTCTCT  
GCTGAACTTGCCTGTTCTGACACCCACATTAATGAGCTGATGATCATTACCATGGGCTTGGTGTCTTA  
GCTGTTCTCTTATGTTGATTGTCTTCTCCTATGTCTGCATTTCTGGGCTGTGTTGGGTATCCCTTCTT  
CTGGAGGGCGATGAAAGCCTTTTCTACCTGTGGCTCCCATCTTACAGTTGTTCTTCTCTATGGATC  
TCTGATGGGTGTGATTTACTTCCCCATCAACTCACTCTACAGAGAGGAAAGTAGGGCTGCTGTTCTC  
TACATGATTGTAATCCCATGCTGAACCCATTCATCTACAGCTTGAGGAACAGAGACATGAAAGAAGCTA  
TGGGTAACCTTTTGGTGGTGGAAAAACAGTTTTCTTATTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR213665 representing NM\_146939  
Red=Cloning site Green=Tags(s)

MGKISRVNQSVASDFLLGLSEQPGEQPLLFGIFLGMVLVTMVGNLLIIFVISSDAHLHTPMYFFLANLS  
 LTDACFTSASVPKMLANIYTSQTIISYSGCLTQLYFLLMFGGLDNCLLAVMAYDRVVAICQPLHYSTAMS  
 PQLCALMLCTCWLTNCPALMHTLLLTRVAFCAHTAIPHFYCDPSALLKLACSDTHINELMIITMGLVFL  
 AVPLMLIVFSYVCISWAVLGISSGGRWKAFTCGSHLTVVLLFYGSLMGVYLLPPSTHSTERESRAAVL  
 YMIVIPMLNPFYISLRNRDMKEAMGKLFGGKTVFLL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9021\\_e07.zip](https://cdn.origene.com/chromatograms/mm9021_e07.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_146939

**ORF Size:** 951 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_146939.1](#), [NP\\_667150.1](#)

**RefSeq Size:** 954 bp

**RefSeq ORF:** 954 bp


**Locus ID:** 258941

**Cytogenetics:** 2 B

**MW:** 35.4 kDa

**Gene Summary:** Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

### Product images:



Circular map for MR213665