

## Product datasheet for **MR213665A1V**

### Mouse Olfr354 (NM\_146939) AAV Particle

#### Product data:

|                           |                                                                             |
|---------------------------|-----------------------------------------------------------------------------|
| Product Type:             | AAV Particles                                                               |
| Product Name:             | Mouse Olfr354 (NM_146939) AAV Particle                                      |
| Tag:                      | Myc-DDK                                                                     |
| Symbol:                   | Olfr354                                                                     |
| Synonyms:                 | MOR127-4                                                                    |
| Mammalian Cell Selection: | None                                                                        |
| Vector:                   | pAAV-AC-Myc-DDK (PS100089)                                                  |
| ORF Nucleotide Sequence:  | >MR213665 representing NM_146939<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGGGAAAAATAAGTAGAGTGAACCAATCTGTTGCTTCAGACTTCCTTCTTCTGGGACTCTCTGAACAGC  
CAGGGGAGCAGCCTCTTCTATTTGGCATCTTCTGGGCATGTACCTGGTCACTATGGTGGGAACCTACT  
CATCATCTTTGTATCAGTTCTGATGCACACCTACATACTCCTATGTACTTCTTCTAGCAAATCTCTCC  
TTAACTGATGCCTGTTTCACTTCTGCATCAGTACCCAAATGCTTGCCAACTTTATACCCAGAGTCAAA  
CCATCTCCTATTCTGGGTGTCTTACTCAGTTGATTTTCTCCTGATGTTTGGCGGCTCGACAACCTGCCT  
TCTGGCTGTGATGGCATATGACCGCTATGTAGCCATCTGCCAGCACTCCATTACAGCACAGCCATGAGT  
CCCCAACTTTGTGATTAATGCTGTGTACATGCTGGGTGCTGACCAACTGCTCCTGCCCTTATGCACACAC  
TGCTGCTGACCCGTGTGGCTTTCTGTGCTCACACAGCCATCCCCATTTCTACTGTGATCCAGTGCTCT  
GCTGAACTTGCTGTTCTGACACCCACATTAATGAGCTGATGATCATTACCATGGGCTTGGTGTCTCTA  
GCTGTTCTCTTATGTTGATTGTCTTCTCCTATGTCTGCATTTCTGGGCTGTGTTGGGTATCCCTTCTT  
CTGGAGGGCGATGGAAGCCTTTTCTACCTGTGGCTCCCATCTTACAGTTGTTCTTCTCTCTATGGATC  
TCTGATGGGTGTGATTTTACTTCCCCATCACTCACTCTACAGAGAGGAAAGTAGGGCTGCTGTTCTC  
TACATGATTGTAATCCCATGCTGAACCATTCATCTACAGCTTGAGGAACAGAGACATGAAAGAAGCTA  
TGGGTAACCTTTTGGTGGTGGAAAAACAGTTTCTTATTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Protein Sequence:</b> | <p>&gt;MR213665 representing NM_146939</p> <p>Red=Cloning site Green=Tags(s)</p> <p>MGKISRVNQSVASDFLLGLSEQPGEQPLLFGIFLGMVLVTMVGNLLIIFVISSDAHLHTPMYFFLANLS<br/> LTDACFTSASVPKMLANIYTSQSTISYSGCLTQLYFLLMFGGLDNCLLAVMAYDRYVAICQPLHYSTAMS<br/> PQLCALMLCTCWLTNCPALMHTLLLTRVAFCAHTAIPHFYCDPSALLKLACSDTHINELMIITMGLVFL<br/> AVPLMLIVFSYVCISWAVLGIPSSGGRWKAFTCGSHLTVVLLFYGSLMGVYLLPPSTHSTERESRAAVL<br/> YMIVIPMLNPFIYSLRNRDMKEAMGKLFGGGKTVFLL</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p> |
| <b>Species:</b>          | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Serotype:</b>         | AAV-2                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>ACCN:</b>             | NM_146939                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>ORF Size:</b>         | 951 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Buffer:</b>           | PBS with 0.001% Pluronic F68                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Stability:</b>        | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.       |
| <b>RefSeq:</b>           | <u>NM_146939.1, NP_667150.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>RefSeq Size:</b>      | 954 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>RefSeq ORF:</b>       | 954 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Locus ID:</b>         | 258941                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Cytogenetics:</b>     | 2 B                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>MW:</b>               | 35.4 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |