

Product datasheet for **MR213699L4V**

Vstm4 (NM_178791) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Vstm4 (NM_178791) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Vstm4
Synonyms:	E130203B14Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_178791
ORF Size:	957 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR213699).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178791.4
RefSeq Size:	2743 bp
RefSeq ORF:	960 bp
Locus ID:	320736
UniProt ID:	T1NXB5
Cytogenetics:	14 B
Gene Summary:	Peptide Lv enhances L-type voltage-gated calcium channel (L-VGCC) currents in retinal photoreceptors.[UniProtKB/Swiss-Prot Function]



[View online »](#)