

Product datasheet for MR213712

Tigit (NM_001146325) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tigit (NM_001146325) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tigit
Synonyms:	ENSMUSG00000071552; Vstm3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR213712 representing NM_001146325 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCATGGCTGGCTGCTCCTGGTCTGGGTCCAGGGGCTGATACAGGCTGCCTTCCTCGCTACAGGAGCCA
 CAGCAGGCACGATAGATACAAAGAGGAACATCTCTGCAGAGGAAGGTGGCTCTGTCATCTTACAGTGTCA
 CTTCTCCTCTGACACAGCTGAAGTGACCAAGTCGACTGGAAGCAGCAGGACCAGCTTCTGGCCATTTAT
 AGTGTGACCTGGGTGGCATGTCGCTTCAGTCTTCAGTGATCGGGTGGTCCAGGCCCCAGCCTAGGCC
 TCACCTTCCAGTCTCTGACAATGAATGACACGGGAGAGTACTTCTGTACCTATCATACGTATCCTGGTGG
 GATTTACAAGGGGAGAATATTCTGAAGGTCCAAGAAAGCTCAGTGGCTCAGTTCCAGACTGCCCCGCTT
 GGAGGAACCATGGCTGCTGTGCTGGGACTCATTTGCTTAATGGTCACAGGAGTGACTGTACTGGCTAGAA
 AGAAGTCTATTAGAATGCATTCTATAGAAAGTGGCCTTGGGAGAACAGAGCGGAGCCACAGGAATGGAA
 CCTGAGGAGTCTCTATCCCCTGGAAGCCCTGTCCAGACACAACTGCCCTGCTGGTCCCTGTGGAGAG
 CAGGCAGAAGATGACTATGCTGACCCACAGGAATACTTTAATGTCCTGAGCTACAGAAGCCTAGAGAGCT
 TCATTGCTGTATCGAAGACTGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR213712 representing NM_001146325
 Red=Cloning site Green=Tags(s)


MHGWLVLVWVQGLIQAFLATGATAGTIDTKRNI SAEEGSVILQCHFSSDTAEVTQVDWKQDQLLAIY
 SVDLGWHVASVFSRVPVGPGLTFQSLTMNDTGEYFCTYHTYPGGIYKGRIFLKVQESSVAQFQTAPL
 GGTMAAVLGLICLMVTGTVLARKKSIRMHSIESGLGRTEAEPQEWNLRLSSPGSPVQTQTAPAGPCGE
 QAEDDYADPQEYFNVLSYRSLESFIAVSKTG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja2276_h02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001146325

ORF Size: 723 bp


OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001146325.1</u> , <u>NP_001139797.1</u>
RefSeq Size:	963 bp
RefSeq ORF:	726 bp
Locus ID:	100043314
Cytogenetics:	16 B4
MW:	26.1 kDa

Product images:



Circular map for MR213712