

## Product datasheet for **MR213811**

### **Olfr1218 (NM\_146818) Mouse Tagged ORF Clone**

#### Product data:

|                           |                                                                                           |
|---------------------------|-------------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                       |
| Product Name:             | Olfr1218 (NM_146818) Mouse Tagged ORF Clone                                               |
| Tag:                      | Myc-DDK                                                                                   |
| Symbol:                   | Olfr1218                                                                                  |
| Synonyms:                 | MOR233-12                                                                                 |
| Mammalian Cell Selection: | Neomycin                                                                                  |
| Vector:                   | pCMV6-Entry (PS100001)                                                                    |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                                      |
| ORF Nucleotide Sequence:  | <p>&gt;MR213811 representing NM_146818</p> <p>Red=Cloning site Blue=ORF Green=Tags(s)</p> |

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGC**

ATGCAAAACCAGAGTTTTGTTACTGAGTTCATAATCTTGGGGCTTCTCAGAACCAATAGTTGAGAAAA  
 TACTGTTTTTGTACTTTTGTAGTCTATCTTGCAACTATAGGGGGCAACATAATAATTGTTGTGACCAT  
 CATGTACAGCCCTGCACTGCTGGGATCCCCATGTATTTCTTCTGGCGTTCTTGCTCTCCTGGATCTA  
 TGTGTCTCTTCCACTGTCATACCAAGATGATTGTAGACTTCTTCTATGAGAAGAAGACCATCTCCTTTG  
 GCTGTTGCATGATGCAACTGTTTTCACTCCATTTCTTCTCAGGAACAGAAATGATTGCTCCTAGCAGCCAT  
 GGCTATGACCGCTACGTGGCTATTTGCAAGCCCTTGCACTATTTTCCATCCTTACCGGAGGCTCTGT  
 AGCATTTTGGTAGCAATATCCTGGGCAGGGGGCTTCTTGCACTCTATCATACAAGTTATATTCACATTGC  
 AGCTGCCTTTATGTGGACCAATGTCTTGTATGATCATTACACATGTGACTTGTTCATTGCTGAAGCTGGC  
 CTGCACGGACACTCATATATTTGTTCTTTTGGTATTTGCAAACAGCGGTGCTATCTGCATCATTATCTTC  
 TCCTTGTTGCTTGTTCCTATGGTGTCAATTTGTTCTCTCTGAGAGCCACAGTTCTGAAGGACGAAGGA  
 AAGCTCTTCCACCTGTGGAGCCACGTTACTGTTGTAGTTTTGTTCTTGTTCCTGTTCCATGCATTAATATA  
 TGCACGGGATACATCTGCATTTTCCTTTGAGAAACACACACTTTTATTTGTCAATGCTCCTGACACCATG  
 CTAATCTACAGTTTACACTTTCAGGAATAAGGAAATGATAAATGCCATCAGAAAAATGTCAAGAGAC  
 TGAATGATTTTTGTTAGATT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

**Protein Sequence:** >MR213811 representing NM\_146818  
 Red=Cloning site Green=Tags(s)


MQNQSFVTEFIILGLSQNPIVEKILFFVLLLVLATIGGNIIVVTIMYSPALLGSPMYFFLAFLSFLDL  
 CVSSTVIPKMI VDFYIEKKTISFGCCMMQLFSVHFFSGTEMIVLAAMAYDRYVAICKPLHYFSILTRRLC  
 SILVAISWAGGFLHSIIQVIFTLQLPLCGPNVIDHYTCDLFPLKLACTDTHIFVLLVFANSGAICIIIF  
 SLLLVS YGVILFSLRAHSSEGRRKALSTCGAHVTVVVFLVPCILIIYARDTSAFSFEKHTLLFVNVLTPL  
 LNPTVYTFRNKEMINAI RKMCKRLKMIFVRF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9069\\_a12.zip](https://cdn.origene.com/chromatograms/mm9069_a12.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_146818

**ORF Size:** 933 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_146818.2](#), [NP\\_667029.2](#)

**RefSeq Size:** 936 bp

**RefSeq ORF:** 936 bp


**Locus ID:** 258815

**Cytogenetics:** 2 E1

**MW:** 35.1 kDa

**Gene Summary:** Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

## Product images:



Circular map for MR213811