

Product datasheet for **MR213834A1V**

Mouse Pnp2 (NM_001123371) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Pnp2 (NM_001123371) AAV Particle
Tag:	Myc-DDK
Symbol:	Pnp2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR213834 representing NM_001123371 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGCAGAAAGGGCAAGATTTGGGCGCCTCTGTCGGTGACACCATAGAGAAAGAGTTCACATACGAAG
ATTATGAGACCACTGCCAAGTGGCTTCTGCAACACACTGAATATCGACCTCAAGTGCCAGTGATCTGTGG
TTCCGGCTTAGGAGGGCTGACTGCTCACTTAAAGGCGGCTCAGATCTTTGACTACAATGAGATACCCAAC
TTTCCCAAAGCACAGTGCAAGGTCACGCAGGCCGACTGGTGTGGATTGCTGAATGGCAGATCCTGTG
TGATGATGCAAGGCCGTTCCATATGTATGAAGGATACTCACTGTCAGAGGTGACATTCCAGTGAGAGT
TTTCCATCTTCTGGGTGTGGAACCTTTGGTGGTCACCAATGCTGCTGGAGGACTCAACCCAATTTTGAA
GTTGGAGATATTATGCTGATCCGTGATCACATCAACCTACCTGGTTTCTGTGGCCAGAACCCTCTCCGGG
GCCCCAACGATGAAAGGTTTGGAGTTCGTTTCTGCCATGTCTGATGCTTATGACCGGGATATGAGGCA
GAAGGCTTTCAGTGCCTGGAAACAAATGGGGGAGCAACGAAAGCTACAAGAAGGCACCTATGTGATGTTG
GCAGGTCCCACTTTGAGACTGTGGCAGAGAGTCGTCTGCTAAAGATGCTGGGGGCAGATGCTGTTGGCA
TGAGCACAGTCCCAGAAGTTATCGTCGCAAGGCACTGTGGGCTTCGTGTCTTTGGTTTCTCACTCATTAC
GAACATGGTTGTCATGGATTATGAGAACTTGAGAAGGCCAATCACAAGGAAGTTCTGGATGCCGGGAAA
GCAGCTGCACAGAAATTGGAAGGTTTGTCTCCATTCTTATGGAGAGCATTCCACTCCCGGAGTGAACC
TCCTT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR213834 representing NM_001123371 Red=Cloning site Green=Tags(s) MEQKGQDLGASVGDITIEKFTYEDYETTAKWLLQHTYRQVAVICGSGLGGLTAHLKAAQIFDYNEIPN FPQSTVQGHAGRLVFGLLNGRSCVMMQGRFHMVEGYSLSEVTFPVRVFHLLGVETLVVTNAAGGLNPNFE VGDIMLIRDHINLPGFCGQNPLRGPNDERFGVRFPMASDAYDRDMRQKAFSAWKQMGQRKLQEGTYVML AGPNFETVAESRLLKMLGADAVGMSTVPEVIVARHCGLRVFGFSLITNMVMDYENLEKANHKEVLDAGK AAQKLERFVSILMESIPLPECNLL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001123371
ORF Size:	915 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001123371.2, NP_001116843.1</u>
RefSeq Size:	1223 bp
RefSeq ORF:	918 bp
Locus ID:	667034
Cytogenetics:	14 C1
MW:	33.9 kDa