

Product datasheet for MR213887

Olfr381 (NM_147022) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr381 (NM_147022) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr381
Synonyms:	MOR135-4; RP23-103M20.1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR213887 representing NM_147022 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGATAATGAATAACAAACTGTCATCACCCAGTTTATTCTCCTGGGCTGCCATCCTGCCAGAGTACC
AGCACCTGTTCTATGCCCTGTTCTGTCCATGTACCTCACCACCATCCTGGGAACCTCATCATCAT
CCTATTCAACTGGACTCCCATCTCCACACCCCATGTACTTGTTCTCAGCAACTTGCCTTCTCTGAT
CTCTGCTTTTCTCTGTGACAATGCCAAGTTGCTGCTGAACATGCAGAGCCAGGACACATCCATCCCT
ATGCAGTTGCTGACACAAATGACTTTTCCAATCTTTTGAAGCCTGGAAATCTTCTCCTTGTGAT
CATGGCCTATGACCGCTATGCGGCCATCTGCCTTCCCCTTCACTACACCAGCATCATGAGCCCCAAGCTC
TGTGTTTGTCTGGTGCTGCTGTCTGGTGATTAGTATGCTGTATTCCATGTTGCACACTTACTCTTGG
CTAGATTGTCATTCTGTGAGGACAATGTGATCCCTCACTTCTTCTGTGACATATCTGCACTGCTCAAGTT
GGCCTGCTCTGACATTATTAATGAACAAATGATATTTTCTTAGGAGGGCCCTTACTGTATCCCCA
TTCTTACTCATTGTTGTGTCCTATATAAAATTGCTTCTCTATTCTAAAGATTTCTACTCGGGCCA
TACACAAGGTCTTCTCCACCTGTGGCTCACACCTATCTGTGGTCTCACTGTTCTATGGGACAATTATTGG
TCTCTACTTATGTCCATCAGCTAATAACTTTTCAGTGAAGAAGGCTTCCATTACCATGATGTACACAGTG
GTGACTCCCATGCTGAACCTTTTCATCTACAGCCTGAGGAACAGAGATATAAAAGAGGCCCTTGTAAAGAG
TCCTTATAAGAAAATATCATTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR213887 representing NM_147022
 Red=Cloning site Green=Tags(s)


MIMNNKTVITQFILLGLPILPEYQHLFYALFLSMYLTITLGNLIIIIQLDShLHTPMYFLSNLSFSD
 LCFSSVTMPKLLLNMQSQDTSIPYAGCLTQMYFSNLFGSLEIFLLVIMAYDRYAAICLPLHYTSIMSPKL
 CVCLVLLSWVISMLYMLHTLLARLSFCEDNVIPHFFCDISALLKLACSDIHINELMIFFLGGPLTVIP
 FLLIVVSYIQIVFSILKISSTRAIHKVFSTCGSHLSVVSIFYGTIIGLYLCPANNFVSVKASITMMYTV
 VTPMLNPFYISLRNRDIKEALVRVLIKKISL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9070_b01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_147022

ORF Size: 933 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_147022.2](#), [NP_667233.2](#)

RefSeq Size: 936 bp

RefSeq ORF: 936 bp


Locus ID: 259024

Cytogenetics: 11 B4

MW: 35.2 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR213887