

Product datasheet for **MR213887A1V**

Mouse Olfr381 (NM_147022) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr381 (NM_147022) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr381
Synonyms:	MOR135-4; RP23-103M20.1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR213887 representing NM_147022 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGATAATGAATAACAAAAGTGCATCACCCAGTTTATTCTCCTGGGCCTGCCATCCTGCCAGAGTACC
AGCACCTGTTCTATGCCCTGTTCTGTCCATGTACCTCACCACCATCCTGGGGAACCTCATCATCATCAT
CCTCATTCAACTGGACTCCCATCTCCACACCCCATGTAATGTTTCTCAGCAACTTGTCTTCTCTGAT
CTCTGCTTTTCTCTGTCACAATGCCAAGTTGCTGCTGAACATGCAGAGCCAGGACACATCCATCCCCT
ATGCAGGTTGCCGACACAAATGTACTTTTCCAATCTTTTGAAGCCTGGAAATCTTCCTCCTTGAT
CATGGCCTATGACCGCTATGCGGCCATCTGCCTTCCCCTTCACTACACCAGCATCATGAGCCCCAAGCTC
TGTGTTTGTCTGGTGTCTGTCTGGGTGATTAGTATGCTGTATTCCATGTTGCACACTCTACTCTTGG
CTAGATTGTCTTCTGTGAGGACAATGTGATCCCTCACTTCTTCTGTGACATATCTGCACTGCTCAAGTT
GGCCTGCTCTGACATTATTAATGAACAAATGATATTTTCTTAGGAGGGCCCTTACTGTATCCCA
TTCTTACTCATTGTTGTCTATATACAAATGTCTTCTCTATTCTAAAGATTTCTACTCGGGCCA
TACACAAGGTCTTCTCCACCTGTGGCTCACACCTATCTGTGGTCTCACTGTTCTATGGGACAATTATTGG
TCTCTACTTATGTCCATCAGTAATAAATTTTCACTGAAGAAGGCTTCCATTACCATGATGTACACAGTG
GTGACTCCCATGCTGAACCTTTCTATCTACAGCCTGAGGAACAGAGATATAAAGAGGCCCTTGAAGAG
TCCTTATAAAGAAAATATCATTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR213887 representing NM_147022 Red=Cloning site Green=Tags(s) MIMNNKTVITQFILLGLPILPEYQHLFYALFLSMYLTITLGNLIIIIILQLDSHLHTPMYFLSNLSFSD LCFSSVTMPKLLLNMQSDTSIPYAGCLTQMYFSNLFGSLEIFLLVIMAYDRYAAICLPLHYTSIMSPKL CVCLVLLSWVISMLYSMLHTLLARLSFCEDNVIPHFCDISALLKLACSDIHINELMIFFLGGPLTVIP FLLIVVSYIQIVFSILKISSTRAIHKVFSTCGSHLSVVSFLFYGTIIGLYLCPSSANNFSVKKASITMMYTV VTPMLNPFIIYSLRNRDIKEALVRVLIKKISL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_147022
ORF Size:	933 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_147022.2, NP_667233.2</u>
RefSeq Size:	936 bp
RefSeq ORF:	936 bp
Locus ID:	259024
Cytogenetics:	11 B4
MW:	35.2 kDa