

## Product datasheet for **MR213890L3V**

### 2310057J18Rik (NM\_026336) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	2310057J18Rik (NM_026336) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	2310057J18Rik
Synonyms:	2210414P04Rik; Leg1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026336
ORF Size:	1011 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR213890).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026336.3</a> , <a href="#">NP_080612.1</a>
RefSeq Size:	1319 bp
RefSeq ORF:	1014 bp
Locus ID:	67719
UniProt ID:	<a href="#">Q8C6C9</a>
Cytogenetics:	10 A4
Gene Summary:	May be involved in early liver development.[UniProtKB/Swiss-Prot Function]



[View online »](#)