

Product datasheet for **MR213893A1V**

Mouse Olfr599 (NM_146731) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr599 (NM_146731) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr599
Synonyms:	MOR23-1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR213893 representing NM_146731 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGCTATACCAACTTGTCATACCTGAACCCAGGGACAGTGATTCTGATTGGCATCCCTGGACTGGAGC
ATGTGCAGTTTTGGATTGGATTTCCCTTCTTTGTTGTATGCCTTGTGGCTTTGCTGGGAAACCTCTTCT
GCTAATCATTATCCCCACAGAACGCAGTCTACACCAACCCATGTACATTTCTTGGCAGTACTCGCAGCT
ACTGACCTTGGTCTCTGCCTAGCCATTGCCCTAAGATGTTGGCCATCTTCTGGTTTGGCTCTGTTCCA
TGGCTTTTGTATGCTTGCCTCACCAGCTCTTCTCATCCATGCCTTGCAAGGAATGGAATCTGGCGTTCT
GCTTGCCATGGCCTTTGACCGATATGTCGCCATCTGCGATCCTCTGAGGCACACAGCCGCTTACACCT
TTATTTTACTTAGAGTCGTTCTGGTGGTTGCAATTCGGGCAACTGTACTCGTTGGGTTTTACCCATTC
TCCTCAAGCGACTGCAATGGTTCCACTCTGTTGTGATTGTTTATTCTTACTGTGAGCACATGGCTGTAGT
CAAGCTGGCGGCAGAAGATGTCCGTATAAATAAATCATATGGGCTCTTTGTGGCTTTGCAATTCTAGGT
TTTGACATGATATTTGTTTTCATCTCCTACATTTGATTTTCCGAGCTGTTTTCTGCTTCTCCTCAGAAAG
AGGCACGAAGCAAAGCATTCAACACTTGTACGGCTCATATTGTTGTCTTCTGGAGTTTTATATCCTTGC
CTTTTTTCTTTTTCAGTCACCGTTTTGGTCACGTATACCGTATGTCCACATCTTATTGTCTACCATC
TATCTACTTCTGCCCCCGCCCTTAACCCATTGTGTATGGTGTGAAGACCAAGGAGATCCGCAATGGG
TTGTCCAAATTTTGTCTGAAGTCTAATACTCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR213893 representing NM_146731 Red=Cloning site Green=Tags(s) MGYTNLSYLNPGTVILIGIPGLEHVQFWIGPFFVCLVALLGNLFLLIIPTERSLHQPMYIFLAVLAA TDLGLCLAIAPKMLAIFWFGSCSMAFDACLTLFFIHALQGMESGVLLAMAFDRYVAICDPLRHTAVLTP LFLLRVVLVVAIRATVLVGVLPIILLKRLQWFHSVVIVHSYCEHMAVVKLAAEDVRINKSYGLFVAFAILG FDMIFVFISYILIFRAVFRLPQKEARSKAFNTCTAHIVVFLEFYILAFFSFFSHRFGHVSPYVHILLSTI YLLLPALNPVYGVKTEIRKWWVQIFVLKSNTQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_146731
ORF Size:	945 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_146731.1, NP_666942.1</u>
RefSeq Size:	948 bp
RefSeq ORF:	948 bp
Locus ID:	258726
Cytogenetics:	7 E3
MW:	35.6 kDa