

Product datasheet for **MR213930A1V**

Mouse Ssxb3 (NM_198898) AAV Particle

Product data:

Product Type: AAV Particles
 Product Name: Mouse Ssxb3 (NM_198898) AAV Particle
 Tag: Myc-DDK
 Symbol: Ssxb3
 Mammalian Cell Selection: None
 Vector: pAAV-AC-Myc-DDK (PS100089)
 ORF Nucleotide Sequence: >MR213930 representing NM_198898
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGAACAGTGAGTTCCTGTGAAAAGATCCCATGGAAGTGCTTTATGAACCAAAGAATATCTGCAAGG
 CCTTCCAGGATATTTCTACATACTTCTCTGATGAAGAGTGGGAAAGCTGACTCAATGGCAGAAAAGTGC
 CTACATGTACATGAAAAGAACTACATCAGAATGACTGACCTAGGGGTACCGTGAACCAACCAGTTTTCT
 ATGCGTGGCAAGGAGCAGGAAAAGCAATCCCTGGTTGAAGGCATTGAAGTTCATGACAGTGATGATGAAT
 GCTTTGAAGGATCTTTTGTGTGACACCGAGAAAACGAATGAAGCTGACGTGAGTACAATAAGTATTCA
 TAATGTAGAAGGAAGCCTTGCCCTGGAAAGAATGATTCCAACCTGGCAGAAAAGTATGGCATTACAGGTC
 AATGTTTGGAGCCACCGCTGCGAGAGAGGAAGTATCGAGTGATATTCAGAGATTAGTGATCTCTGAAG
 AGGAAGAAGATGATGATTAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR213930 representing NM_198898
 Red=Cloning site Green=Tags(s)

METVSSCEKIPMEVL YEPKNICKAFQDISTYFSDEEWGKLTQWQKSAYMYMKRNYIRMTDLGVTVNQPVF
 MRGKEQEKQSL VEGIEVHDSDECFEGSFDVTPRKRMKLT SVTISIHNVEGSLASGKNDSNLAETDGIQV
 NVWSHRLRERKYRVIYSEISDP EEEEEDDY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
 Serotype: AAV-2
 ACCN: NM_198898



[View online »](#)

ORF Size:	510 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_198898.2</u>
RefSeq Size:	604 bp
RefSeq ORF:	513 bp
Locus ID:	278174
Cytogenetics:	X A1.1
MW:	19.8 kDa