

Product datasheet for **MR213947A1V**

Mouse Olfr871 (NM_146903) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr871 (NM_146903) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr871
Synonyms:	MOR141-2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR213947 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAACCAGAAAACCAACAAATATTTTAGAATTTTCTTCTGGGGTTTTCCCAATATCCAGAGCATC
AACCCATGCTATTTGGACTGTTTCTGCTCATGTTTGGTTCGCTGTGTTGGGAATCTTCTCATCATTCT
GGCCGTCAGCATTGACTCTCACCTGCATACTCCCATGTACTTCTTCTATCTAACCTGTCCTTTTCTGAC
ATTGGTTTCATCTCTACAACATCCCTAAGATGTTGGTGAATATCCAAACACAGAGCAAGTCCATCTCCT
ATGCAGAATGCATCACCAGATTTATTTTTCATGCTCTTTGGAGGCATGGACATACTTCTCCTCACCGT
GATGGCCTATGACCGATTTGTGGCCATCTGTCACCCCCTTCACTATTCACTGTCATTATGAATCCCCAACTA
AGTGGCTTGCTGGTTCTTGATCATGGTTTATTAGCTTTTCATATTCTCTGATACAGAGTCTATTGATGC
TGCGGTTGTCCTTCTGTACAAATCAGATAATTAACACTTTTACTGTGAATATTCTAGAGCCCTCACTAT
AGCCTGCTCAGACACACTAATCAATCATATCCTTCTTTATATTCTGATATGTGTCCTTGGCTTCATCCCT
TTCTCAGGGATCCTTTATTCACTGTAAAATTGTTTCTTCTATTTTGAGAATTCATCAACAGATGGAA
AATATAAAGCATTTTCTACCTGTGGTCTCATCTATCAGTGGTTTCTTTATTCTATGGGACAGGCCTTGG
TGTGTACCTTAGTTCTGATGTAACCTCCTCTGGGAAGGACGTGGTGCCTCAGTAATGTATACAGTG
GTCACCCCTATGCTGAACCTTTTCATCTACAGCTTGAGGAACAAAGACATTAAGAAGGCCTTAAAAAC
TTGGGAGAATACTTCTTTAAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR213947 protein sequence Red=Cloning site Green=Tags(s) MKPENQTNILEFLLLGFSDQPEHQPMFLFGLFLLMFVAVFGNLLIILAVSIDSHLHTPMYFFLSNLSFSD IGFISTTIPKMLVNIQTQSKSISYAECITQIYFFMLFGGMDILLTVMAYDRFVAICHPLHYSVIMNPQL SGLLVLSWFIISYSYSLIQSLLMLRSLFCTNQIIKHFYCEYSRALTACSDTLINHILLYILICVLGFIP FSGILYSYCKIVSSILRIPSTDGKYKAFSTCGSHLSVVSFLFYGTGLGVYLSSDVTSSSGKDVVASVMYTV VTPMLNPFIIYSLRNKDIKKALKTLGRILLK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_146903
ORF Size:	933 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_146903.2</u> , <u>NP_667114.2</u>
RefSeq Size:	1147 bp
RefSeq ORF:	936 bp
Locus ID:	258905
Cytogenetics:	9 A3
MW:	35 kDa