

Product datasheet for MR213965L3V

Slc38a6 (NM_001037717) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc38a6 (NM_001037717) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc38a6
Synonyms:	AW322671; EG625098
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001037717
ORF Size:	1371 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR213965).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001037717.3 , NP_001032806.2
RefSeq Size:	2959 bp
RefSeq ORF:	1374 bp
Locus ID:	625098
UniProt ID:	G3UWV3
Cytogenetics:	12 C3
Gene Summary:	Probable sodium-dependent amino acid/proton antiporter, could be a neuronal transporter for glutamate.[UniProtKB/Swiss-Prot Function]



[View online »](#)