

## Product datasheet for MR213972

### Krtap11-1 (NM\_001113406) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Krtap11-1 (NM\_001113406) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Krtap11-1  
**Synonyms:** AA162070; AA589467; AV018433; Hacl1; KAP11.1; PRED29  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR213972 representing NM\_001113406  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGTCTTTCAACTGCTCCACAAGAAATTGCTCTTCCAGGCCAGTTGGAGGACGTTAACTGCTCCAGTGG  
 GCCAGTTACCACAGCCTCAGCCCGTGACGCAGACTGCCTGAGTGGCCTCTACTTGCCAGTTCTTTCCA  
 AACTGGCTCCTGGCTCCTGGACCATTGTCAGGAGTCATACTGTGAGCCTACTGTCTGCCAGCCAACGTGC  
 TACCAGAGAACTTCTTGATCTCCACTCCTGCCAGGTAACCTGCAATCGACAGACTACCTGCGTCTCCA  
 ATCCTTGCTCTAGCCCTGTAGCCGGCCACTCACTTTTGTGTCCACTGGCTGTCAGCCCTGGGAGGGAT  
 TTCCAGCTCCTGCCAACCAGTGGGAGGCATCTCTACCACCTGCCAACCAGTGGGAGGCATCTCTACCACC  
 TGCCAACAAGTGGGAGGCATCTCTACCGTCTGCCAACCAGTGGGAGGTATCTCTACTGTCTGCCAACCA  
 CCTGTGGAGTCTCCAGGACACACCAGCAGTCTGTGTATCCAGCTGCCAAGAAGTTC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR213972 representing NM\_001113406  
 Red=Cloning site Green=Tags(s)

MSFNCSTRNCSRPVGGRYTAPVGPVTTASARDADCLSGLYLPSSFQTSWLLDHCQESYCEPTVCQPTC  
 YQRTSCISTPAQVTCNRQTTCVSNPCSTPCSRPLTFVSTGCQPLGGISSCQPVGGISTTCQPVGGISTT  
 CQVGGISTVCQPVGGISTVCQPTCGVSRTHQQSCVSSCRRTC

**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001113406

**ORF Size:** 549 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001113406.1](#), [NP\\_001106877.1](#)

**RefSeq Size:** 1008 bp


**RefSeq ORF:** 552 bp

**Locus ID:** 16693

**Cytogenetics:** 16 C3.3

**MW:** 19 kDa

Product images:



Circular map for MR213972